

NEB

expressions

a scientific update



in this issue

- 2** Improving enzymatic DNA fragmentation for next generation sequencing library construction
- 6** NEBNext® Ultra™ II FS DNA Library Prep Kit
- 8** NEBNext® Ultra™ II RNA Library Prep
- 9** Programmable CRISPR/Cas9 Nucleases

Did you get the message?



- 10** Everything you need for your RNA-related workflows:
NEB's murine RNase Inhibitor
- 11** High Yield RNA Synthesis in all Flavors
- 12** Monarch Total RNA Miniprep Kit
- 14** RNA Quantitation by RT-qPCR: Latest enzyme technology for superior RT-qPCR Performance
- 16** NEB's DNA ladders for your daily lab work

Improving Enzymatic DNA Fragmentation for Next Generation Sequencing Library Construction

by Lynne Apone, Ph.D., Eileen Dimalanta, Ph.D., and Fiona Stewart, Ph.D., New England Biolabs, Inc.

INTRODUCTION

The Human Genome Project (HGP), which was officially completed in 2003, was considered to be one of the world's largest collaborative projects of its time (1). This involved many research groups worldwide and had the lofty goal of deciphering all 3 billion bases of the human genome. The project cost almost \$4 billion dollars and took 13 years to complete with the available technology. Over a decade later, advancements in next generation sequencing (NGS) technologies have enabled sequencing of a human genome to become routine, taking less than two days, and at a tiny fraction of the cost of the original HGP.

The ability to quickly and inexpensively sequence whole genomes has truly revolutionized genomics research. Where once single genes or families of genes were studied, now whole genomes, exomes, transcriptomes and epigenomes are interrogated. With recent advances, such as the ability to multiplex and sequence many samples at once, NGS has transitioned from a basic research tool into the clinic, where it impacts discovery, diagnostics and treatment of disease.

Advances in genomics driven by NGS, as well as advances in the technology itself, continue at an amazing pace and move us closer to the realization of personalized medicine, where clinical decisions are tailored to an individual's genome. However, if this pace is to continue, advances in all aspects of the technology must also continue. This includes early steps of the sequencing workflow, specifically in the preparation of samples, before they are sequenced.

To date, there are no sequencing platforms that can sequence intact DNA. Therefore, prior to sequencing, DNA molecules must be fragmented, or broken, into smaller pieces. These DNA fragments are then converted into libraries, by different methods depending on the sequencing platform to be used (Figure 1). In all cases, the libraries generated consist of the fragments of the unknown DNA to be sequenced, flanked by pieces of known DNA (adaptors), which are specific to each sequencing platform.

DNA FRAGMENTATION APPROACHES

One of the major bottlenecks to sample prep is the first step: DNA fragmentation.

The size of the DNA fragments generated depend on the sequencing platform being used, and can range from several hundred base pairs



FIGURE 1: Traditional library preparation workflow



for short read sequencing technologies (e.g., Illumina®, Ion Torrent™) to >10 kb pieces for long read sequencing technologies (e.g., Pacific Biosciences® and Oxford Nanopore Technologies®). Methods for fragmenting DNA are broadly split into two basic categories: mechanical and enzyme-based. Mechanical shearing methods include acoustic shearing, hydrodynamic shearing and nebulization, while enzyme-based methods include transposons, restriction enzymes and nicking enzymes. Although many different options exist to fragment DNA, final fragment size, amount of starting material, upfront capital investment, and scalability must be considered when choosing a fragmentation method. Critically, in order to be useful for NGS, the method used must shear the DNA sufficiently randomly, so that the library being sequenced is fully representative of the original sample.

Mechanical Shearing

Options for mechanical fragmentation of DNA range from small plastic nebulizer devices to sophisticated electronic instruments. The most commonly used technique utilizes **focused**

acoustic shearing devices, such as the instruments made by Covaris®. This involves focused transmission of high-frequency, short wavelength acoustic energy on the DNA sample. The size of the DNA fragments generated (150 – 5,000 bp) is controlled by changing both the intensity and the duration of the acoustic waves, and the protocols used are the same regardless of the amount or GC content of the DNA. Cost, challenges of scalability and sample loss (often caused by sample transfer after shearing), are some of the reasons that users of this method seek alternatives, especially as throughput increases.

If larger DNA fragments are required, **hydrodynamic shearing** can be used. In this method, hydrodynamic shear forces are applied by pushing DNA through the small orifice of a syringe. Size is controlled by altering the speed at which the DNA is pushed through the syringe. Centrifugation can also be used to create hydrodynamic force, by pulling the DNA sample through a hole with a defined size. Here, the rate of centrifugation determines the degree of DNA fragmentation. DNA fragments generated with hydrodynamic shear forces are typically in the

DNA fragmentation landscape prior to NEBNext Ultra II FS

Focused acoustic shearing

- +**
 - Historical "gold standard"
 - Simple control of fragment size
 - Same protocol regardless of GC content or input amount
- - Significant upfront instrument cost
 - Significant disposables cost and waste generation (tubes/plates)
 - Time-consuming
 - Challenging to scale
 - Inflexible reaction volumes
 - Prone to DNA damage

Enzymatic fragmentation

- +**
 - Scalable
 - Flexible reaction volumes
 - Ease of use
 - Minimized DNA damage
- - Prone to sequence bias
 - May be sensitive to contaminants in the DNA sample
 - Significant optimization may be required
 - Transposases can have very specific input requirements



range of 1-75 kb, but require large DNA input amounts (> 1 µg) and throughput is low.

Nebulization is another method used to mechanically fragment DNA. Nebulization uses compressed air to force DNA through a small hole in a nebulizer unit and DNA fragment size is determined by the pressure used. Although this method is inexpensive and fragment size is somewhat tunable (typically 700 – 5000 bp in size), microgram quantities of DNA are required for starting material, and the method is most suitable for small numbers of samples.

Enzymatic Fragmentation

Enzyme-based fragmentation of DNA is an attractive alternative to mechanical shearing methods, as it does not require upfront capital expense, is amenable to quickly processing many samples at the same time, and reduces sample loss. Historically, the main concern with this method has been sequence bias, as many enzymes that cleave DNA have recognition sequences or sequence preferences.

Transposases fragment DNA by cleaving and inserting a short double-stranded oligonucleotide to the ends of the newly cleaved molecule. The inserted oligonucleotide must contain a sequence that is specific to the particular transposase being used. While this method is fast and has low input requirements, the known sequence bias associated with transposases make them incompatible with some applications.

The great majority of **restriction enzymes** have very specific recognition and cleavage sites, and therefore are not suitable for the random cleavage required for most NGS applications. However, restriction enzymes do have utility in a workflow called RAD-Seq (Restriction-site Associated DNA Sequencing). Here, the sequence bias from the restriction enzyme cut site is exploited to target certain regions for sequencing.

Lastly, **non-specific nicking enzymes** can be used to fragment DNA. These enzymes have less sequence bias than transposases or restriction enzyme-based methods, and generate fragments of different sizes (generally 50 bp-1 kb) in a time-dependent manner: the longer the reaction time, the smaller the fragments. Historically, these enzymes have required significant reaction optimization. They have been sensitive to the buffering conditions of the DNA sample, and required different reaction conditions for different DNA input amounts and for varying GC content. While appealing in theory, this has made enzymatic fragmentation methods more challenging to implement, especially in laboratories where a variety of sample types and amounts are used.

Emerging drawbacks to the gold standard

Of all the techniques described, acoustic shearing has traditionally been the method of choice for short read sequencing technologies, such as Illumina; its popularity a result of robust shearing

FIGURE 2:
NEBNext Ultra II FS Kit workflow



with low bias. However, this method requires a significant financial investment in equipment, which can be prohibitive for many researchers. In addition, throughput is low and sample loss high, creating bottlenecks and limitations for users of the technology. Interestingly, recent studies have also shown that oxidative damage can occur during mechanical shearing processes (2,3).

IMPROVING DNA FRAGMENTATION FOR NGS LIBRARY CONSTRUCTION

To address the challenges associated with existing fragmentation approaches, NEB has developed a fragmentation system, the NEBNext® Ultra™ II FS DNA Library Prep Kit (NEB #E7805, E6177), in which a unique enzymatic fragmentation reagent is fully integrated into library preparation to generate low bias, high quality NGS libraries, with a simple, scalable workflow (for more information see page 6).

In order to reduce the NGS sample prep bottleneck, improvements in both performance and ease of use were necessary. In this work, we have focused on the DNA fragmentation step. Our new DNA fragmentation reagent is combined with end repair and dA-tailing reagents, and subsequent adaptor ligation is also carried out in the same vial (Figure 2). For low input samples, PCR amplification is performed prior to sequencing.

Importantly, enzymatic shearing of DNA with this method does not introduce bias into the library, and this method is suitable for input DNA with

a full range of GC content. Reduced sample loss and increased efficiencies of the workflow enable use of lower input amounts, with a range of 100 pg – 0.5 µg, and insert sizes of 100 bp to 1 kb can easily be generated.


PERFORMANCE

Increased Library Yields

The use of enzymatic fragmentation can result in higher library yields as compared to mechanical shearing workflows, as the latter results in sample loss and DNA damage. Achieving sufficient library yields for high quality sequencing from very low input amounts can be especially challenging with mechanical shearing of DNA, a situation compounded by the preference to amplify libraries using as few PCR cycles as possible. Integration of our unique fragmentation reagent with end repair and dA-tailing, removing sample cleanup prior to ligation and eliminating multiple transfer steps all help to minimize sample loss. When combined with the high reaction efficiencies of each step in the workflow and lack of DNA damage caused by mechanical shearing, NEBNext Ultra II FS generates higher yields than library preparation using mechanical shearing methods. High library yields can be achieved with input amounts as low as 100 pg of human genomic DNA with amplification, or as low as 50 ng for PCR-free workflows (Figure 3, page 3).

FIGURE 3:
NEBNext Ultra II FS DNA produces higher yields of PCR-free libraries

Libraries were prepared from Human NA19240 genomic DNA using the input amounts shown, without amplification. NEBNext Ultra II FS libraries were prepared using a 20-minute fragmentation time. "Covaris" libraries were prepared by shearing each input amount in 1X TE Buffer to ~200 bp using a Covaris instrument, followed by library construction using the NEBNext Ultra II DNA Library Prep Kit (NEB #E7645). Library yields were determined by qPCR using the NEBNext Library Quant Kit for Illumina (NEB #E7630). Error bars indicate standard deviation for an average of 2 libraries.



Improvements in Library Quality

As mentioned above, a historical concern regarding the use of enzymatic fragmentation methods was the potential for introduction of bias into a sequencing library, and ultimately into the sequencing data. As shown, this new fragmentation reagent provides consistent uniformity of GC coverage at the full range of input amounts (Figure 4) and over a broad range of GC composition (Figure 5).

In addition, oxidative damage markers typically associated with mechanical-shearing methods (2,3) are observed at significantly lower frequency

in libraries made using this new fragmentation system (Figure 6, page 5). Importantly, a greater difference between the observance of these markers is seen with lower input amounts. This highlights the higher quality of libraries constructed with this new method compared to Covaris-sheared DNA libraries, especially at low input amounts.

EASE OF USE

Robustness of DNA Fragmentation:

Consistent and reliable fragmentation is critical for a new method to be adopted. We optimized

this new fragmentation system to be insensitive to variables such as input amounts, GC content, and DNA buffer conditions. In practice, these details are often unknown for a sample, requiring clean up and quantification prior to traditional enzymatic DNA fragmentation methods. Even when all of the variables are known, traditional enzymatic methods require different fragmentation parameters for each type of sample and DNA input amount. This new fragmentation system addresses all of these issues by requiring just a single-fragmentation protocol for the full range of input amounts (100 pg – 0.5 µg) (Figure 7, page 5) and for the full range of GC content (Figure 5). Also, input DNA can be in water, Tris, 0.1X TE or 1X TE (Figure 8, page 5). Fragmentation using the new system is time dependent, and final library sizes ranging from 100 bp – 1 kb can be generated by simply changing the fragmentation time.

CONCLUSION:

The continued expansion of the use of next generation sequencing depends in significant part on overcoming the limitations and bottlenecks associated with high-quality library preparation, including the initial DNA fragmentation step. While acoustic shearing has for some time been the method of choice for NGS, limitations in terms of instrumentation, throughput and sample damage necessitate sourcing an alternative solution for many users. This new method for enzymatic DNA fragmentation and library

FIGURE 4:
NEBNext Ultra II FS DNA provides uniform GC coverage with human DNA over a broad range of input amounts

Libraries were prepared from Human NA19240 genomic DNA using the input amounts and number of PCR cycles shown, and a 20-minute fragmentation time was used. Libraries were sequenced (2 x 76 bp) on an Illumina® MiSeq®. Reads were mapped to the hg19 reference genome using Bowtie 2.2.4 and GC coverage information was calculated using Picard's CollectGCBiasMetrics (v1.117). Expected normalized coverage of 1.0 is indicated by the horizontal grey line, the number of 100 bp regions at each GC% is indicated by the vertical grey bars, and the colored lines represent the normalized coverage for each library.




FIGURE 5:
NEBNext Ultra II FS DNA provides uniform GC coverage for microbial DNA over a broad range of GC composition

Libraries were prepared using 1 ng of a mix of genomic DNA samples from *Haemophilus influenzae*, *Escherichia coli* (K-12 MG1655) and *Rhodospseudomonas palustris*, with 9 PCR cycles, and sequenced on an Illumina MiSeq. NEBNext Ultra II FS libraries were prepared using a 20-minute fragmentation time. "Covaris" libraries were prepared by shearing 1 ng of DNA in 1X TE Buffer to an insert size of ~200 bp using a Covaris instrument, followed by library construction using the NEBNext Ultra II DNA Library Prep Kit (NEB #E7645). Reads were mapped using Bowtie 2.2.4 and GC coverage information was calculated using Picard's CollectGCBiasMetrics (v1.117). Expected normalized coverage of 1.0 is indicated by the horizontal grey line, the number of 100 bp regions at each GC% is indicated by the vertical grey bars, and the colored lines represent the normalized coverage for each library.








FIGURE 6:
NEBNext Ultra II FS DNA libraries show reduced markers of oxidative damage compared to libraries produced by mechanical shearing

Libraries were prepared from 1 ng and 100 ng Human NA19240 genomic DNA, using 9 and 4 PCR cycles, respectively. NEBNext Ultra II FS libraries were prepared using a 20-minute fragmentation time. "Covaris" libraries were prepared by shearing each input amount in 1X TE Buffer to ~200 bp using a Covaris instrument, followed by library construction using the NEBNext Ultra II DNA Library Prep Kit (NEB #E7645). Libraries were sequenced on an Illumina® HiSeq® 2500 (2 x 75 bp). 723M reads were randomly sampled (seqtk) and aligned to the GRCh38 full reference genome using bwa (0.7.15). Adaptors were trimmed prior to alignment using trimadap (r9). Duplicates were marked using sambalster (0.1.24). Variants were called on chromosome 1 using freebayes (1.0.2.29) with frequency based options requiring at least 10 reads per site. More variants are observed for C>A and G>T transversions compared with all other variants for PCR-amplified Covaris libraries. These variants indicative of oxidative damage are not pronounced in NEBNext Ultra II FS libraries.



preparation addresses these issues, further streamlines the process and improves the quality of NGS libraries. The broadening of the input amount range to as low as 100 picograms enables access to high-quality sequencing of samples not achievable by previous methods, and the reliability and ease of use of the method not only allows automation, but also successful adoption by users with a wide range of laboratory skills.

References:

1. <https://www.battelle.org/docs/default-source/misc/battelle-2011-misc-economic-impact-human-genome-project.pdf>
2. M. Costello et al. (2013) *Nucleic Acids Research*, 41, e67.
3. L. Chen, et al. (2017) *Science*, 17, 355: 752-756.



FIGURE 7:
NEBNext Ultra II FS DNA provides consistent fragmentation regardless of input amount

Libraries were prepared from Human NA19240 genomic DNA using the input amounts shown. NEBNext Ultra II FS libraries were prepared using a 20-minute fragmentation time. Library size was assessed using the Agilent® Bioanalyzer®. Low input (1 ng and below) libraries were loaded on the Bioanalyzer without a dilution. High input libraries were loaded with a 1:5 dilution in 0.1X TE.





FIGURE 8:
NEBNext Ultra II FS DNA provides consistent fragmentation of DNA in water, Tris or TE

Libraries were made using 100 ng Human NA19240 genomic DNA using the NEBNext Ultra II FS kit or the Kapa HyperPlus Kit. Fragmentation conditions targeting ~200 bp inserts were used, which would generate ~320 bp libraries. For the NEBNext Ultra II FS kit, input DNA was in H₂O, Tris, 0.1X TE or 1X TE. For the Kapa HyperPlus kit, libraries were made using the recommended dilution of the supplied Conditioning Solution (CS), or using DNA in Tris, 0.1X TE or 1X TE, in the absence of either Conditioning Solution or 3X bead clean up. Library size distribution was assessed using the Agilent Bioanalyzer. Fragmentation is consistent for the NEBNext Ultra II FS kit for DNA in H₂O, Tris, 0.1X TE or 1X TE.



You'll be thrilled to pieces.



Novel Enzymatic Fragmentation System


Do you need a faster or more reliable solution for DNA fragmentation and library construction? Our new NEBNext Ultra II FS DNA Library Prep Kit meets the dual challenges of constructing high quality next gen sequencing libraries from ever-decreasing input amounts, and simple scalability. Further, this kit provides a fragmentation and library prep solution that can handle different input amounts and sample types, quickly and reliably.

The Ultra II FS kit includes a new fragmentation reagent, which is also combined with end repair and dA-tailing reagents, enabling these steps to be performed in the same tube, with no clean-up steps or sample loss. The same fragmentation protocol is used for any input amount (100 pg–500 ng), or GC content.

You'll be thrilled to pieces with the result – a reliable, flexible, high-quality library prep that is fast and scalable.

FIGURE 1:
NEBNext Ultra II FS DNA produces the highest yields, from a range of input amounts

Libraries were prepared from Human NA19240 genomic DNA using the input amounts and numbers of PCR cycles shown. For NEBNext Ultra II FS, a 20-minute fragmentation time was used. For Kapa™ HyperPlus libraries, input DNA was cleaned up with 3X beads prior to library construction, as recommended, and a 20-minute fragmentation time was used. Illumina® recommends 50 ng input for Nextera®, and not an input range; therefore, only 50 ng was used in this experiment. "Covaris®" libraries were prepared by shearing each input amount in 1X TE Buffer to an insert size of ~200 bp using a Covaris instrument, followed by library construction using the NEBNext Ultra II DNA Library Prep Kit (NEB #E7645). Error bars indicate standard deviation for an average of 3–6 replicates performed by 2 independent users.



View more performance data by downloading our technical notes at NEBNextUltraII.com



High-yield, Scalable Library Preparation with the NEBNext Ultra II FS DNA Library Prep Kit



Next generation DNA library construction (LC) for high-throughput genomics – Data presented by Peter Ellis, Senior Staff Scientist at the Wellcome Trust Sanger Institute

Advantages

- Perform fragmentation, end repair and dA-tailing with a **single enzyme mix**
- Experience **reliable fragmentation with a single protocol**, regardless of DNA input amount or GC content
- Prepare high quality libraries from a **wide range of input amounts**: 100 pg–500 ng
- Use with DNA in standard buffers (e.g., TE, Tris-HCl) and water
- Save time with a **streamlined workflow**: ~ 2.5 hours, with < 15 minutes hands-on time
- Experience **reliable fragmentation**, even with very low input amounts
- Generate **high yields** with increased reaction efficiencies and minimized sample loss
- Vary incubation time to generate a **wide range of insert sizes**



For more information and to request a sample, visit NEBNextUltraII.com or contact your local distributor!

NEBNext® Ultra™ II FS DNA Library Prep Kit

NEW PRODUCTS

FIGURE 2:

NEBNext Ultra II FS DNA provides uniform GC coverage for microbial genomic DNA over a broad range of GC composition

Libraries were prepared using 1 ng of a mix of genomic DNA samples from *Haemophilus influenzae*, *Escherichia coli* (K-12 MG1655), *Rhodopseudomonas palustris* and the library prep kits shown with 9 PCR cycles for consistency across samples, and sequenced on an Illumina MiSeq. NEBNext Ultra II FS libraries were prepared using a 20-minute fragmentation time. For Kapa HyperPlus libraries, input DNA was cleaned up with 3X beads prior to library construction, as recommended, followed by a 25-minute fragmentation time. "Covaris" libraries were prepared by shearing 1 ng of DNA in 1X TE Buffer to an insert size of ~200 bp using a Covaris instrument, followed by library construction using the NEBNext Ultra II DNA Library Prep Kit (NEB #E7645). Reads were mapped using Bowtie 2.2.4 and GC coverage information was calculated using Picard's CollectGCBiasMetrics (v1.117). Expected normalized coverage of 1.0 is indicated by the horizontal grey line, the number of 100 bp regions at each GC% is indicated by the vertical grey bars, and the colored lines represent the normalized coverage for each library.




TABLE 1:
NEBNext Ultra II FS DNA Library Prep Workflow

| | Fragmentation/ End Repair/ dA-Tailing | Adaptor Ligation | Clean Up/ Size Selection | Amplification | Clean Up | Total Workflow |
|--------------|---------------------------------------------|---------------------|-----------------------------|---------------|-----------|-------------------|
| Hands-On | 2 min. | 1 min. | 5 min. | 0–1 min. | 0–5 min. | 8–14 min. |
| Total | 37–62 min. | 16–31 min. | 27–37 min. | 0–34 min. | 0–27 min. | 1.3–3.2 hr. |

ORDERING INFORMATION

| PRODUCTS | NEB # | SIZE |
|---------------------------------------------------------------------|----------|------------|
| NEBNext Ultra II FS DNA Library Prep Kit for Illumina | E7805S/L | 24/96 rxns |
| NEBNext Ultra II FS DNA Library Prep with Sample Purification Beads | E6177S/L | 24/96 rxns |
| NEBNext Ultra II FS DNA Module | E7810S/L | 24/96 rxns |

| ALSO AVAILABLE | NEB # | SIZE |
|----------------------------------------------------------------------|-------------------------------|-------------|
| NEBNext Multiplex Oligos for Illumina (Index Primers Set 1, 2, 3, 4) | E7335, E7500, E7710, E7730S/L | 24/96 rxns |
| NEBNext Multiplex Oligos for Illumina (96 Index Primers) | E6609S/L | 96/384 rxns |
| NEBNext Multiplex Oligos for Illumina (Dual Index Primers Set 1) | E7600S | 96 rxns |

Available products:

• NEBNext Ultra II FS DNA Library Prep Kit for Illumina

Includes optimized mixes for DNA library preparation (enzymatic fragmentation/end repair/dA-tailing, adaptor ligation and PCR enrichment steps) for sequencing on the Illumina platform. This kit includes a new DNA fragmentation reagent, which is also combined with end repair and dA-tailing reagents, enabling these steps to be performed in the same tube, with no clean-ups or sample loss.

NEB # E7805S/L **24/96 rxns**

• NEBNext Ultra II FS DNA Library Prep with Sample Purification Beads

Includes optimized mixes for DNA library preparation (enzymatic fragmentation/end repair/dA-tailing, adaptor ligation and PCR enrichment steps) plus SPRIselect® beads for size selection and cleanup. This kit includes a new DNA fragmentation reagent, which is also combined with end repair and dA-tailing reagents, enabling these steps to be performed in the same tube, with no clean-ups or sample loss.

NEB # E6177S/L **24/96 rxns**

• NEBNext Ultra II FS DNA Module

This module is part of the NEBNext Ultra II FS workflow, and includes optimized mixes for DNA library preparation (enzymatic fragmentation/end repair/dA-tailing, adaptor ligation and PCR enrichment steps).

NEB # E7810S/L **24/96 rxns**


Easy, Fast & Extremely Efficient Library Prep Kits for RNA-Seq

NEBNext® Ultra™ II RNA Library Prep

Do you need increased sensitivity and specificity from your RNA-seq experiments? Do you have ever-decreasing amounts of input RNA? To address these challenges, our next generation of RNA library prep kits have been reformulated at each step, resulting in several fold higher yields of high quality libraries, enabling use of lower input amounts and fewer PCR cycles. NEBNext Ultra II RNA Kits have streamlined, automatable workflows and are available for directional (strand-specific, using the “dUTP method”^{1,2}) and non-directional library prep, with the option of SPRISelect® beads for size selection and clean-up steps.

NEBNext Ultra II Directional RNA produces the highest yields, from a range of input amounts

Poly(A)-containing mRNA was isolated from Universal Human Reference RNA (Agilent® #740000) and libraries were made using the NEBNext Ultra II Directional RNA Kit (plus the NEBNext poly(A) mRNA Magnetic Isolation Kit), Kapa Stranded mRNA-Seq Kit, Kapa mRNA HyperPrep Kit and Illumina TruSeq Stranded mRNA Kit. The input RNA amount and number of PCR cycles are indicated. Library yields from an average of three replicates are shown.



1. Parkhomchuck, D., et al. (2009) *Nucleic Acids Res.* 37, e123.
2. Levin, J. Z., et al. (2010) *Nature Methods* 7, 709–715.

ORDERING INFORMATION

| PRODUCT | NEB # | SIZE |
|------------------------------------------------------------------------------|----------|------------|
| NEBNext Ultra II Directional RNA Library Prep Kit for Illumina | E7760S/L | 24/96 rxns |
| NEBNext Ultra II Directional RNA Library Prep with Sample Purification Beads | E7765S/L | 24/96 rxns |
| NEBNext Ultra II RNA Library Prep Kit for Illumina | E7770S/L | 24/96 rxns |
| NEBNext Ultra II RNA Library Prep with Sample Purification Beads | E7775S/L | 24/96 rxns |

OTHER PRODUCTS YOU MIGHT BE INTERESTED IN:

| | | |
|---------------------------------------------------------------------------------|-------------------------------|--------------|
| NEBNext rRNA Depletion Kit (Human/Mouse/Rat) | E6310S/L/X | 6/24/96 rxns |
| NEBNext rRNA Depletion Kit (Human/Mouse/Rat) with RNA Sample Purification Beads | E6350S/L/X | 6/24/96 rxns |
| NEBNext Poly(A) mRNA Magnetic Isolation Module | E7490S/L | 24/96 rxns |
| NEBNext Ultra II RNA First Strand Synthesis Module | E7771S/L | 24/96 rxns |
| NEBNext Ultra II Directional RNA Second Strand Synthesis Module | E7550S/L | 24/96 rxns |
| NEBNext Ultra II Non-Directional RNA Second Strand Synthesis Module | E6111S/L | 20/100 rxns |
| NEBNext Multiplex Oligos for Illumina (Index Primers Set 1, 2, 3, 4) | E7335, E7500, E7710, E7730S/L | 24/96 rxns |
| NEBNext Multiplex Oligos for Illumina (96 Index Primers) | E6609S/L | 96/384 rxns |
| NEBNext Multiplex Oligos for Illumina (Dual Index Primers Set 1) | E7600S | 96 rxns |

Advantages

- **Generate high yield**, high-quality libraries even with limited amounts of RNA:
 - 10 ng – 1 µg total RNA (poly(A) mRNA workflow)
 - 5 ng – 1 µg total RNA (rRNA depletion workflow)
- **Minimize bias**, with fewer PCR cycles required
- **Increase library complexity** and transcript coverage
- **Increase flexibility** by ordering reagents specific to your workflow needs
 - Directional and non-directional kits available
 - rRNA depletion and poly(A) mRNA isolation reagents supplied separately
 - Adaptors and primers (12-, 96-, and dual index) supplied separately
- Enjoy the **reliability of the gold standard SPRISelect size selection and clean-up beads**, supplied in just the amounts you need
- Save time with **streamlined workflows**, reduced hands-on time, and automation compatibility
- Rely on **robust performance**, even with low quality RNA, including FFPE



For more information and to request a sample, visit **UltraIRNA.com** or contact your local distributor!

Gene Editing

NEW PRODUCTS


Programmable CRISPR/Cas9 Nucleases

NEB now offers **EnGen® Spy Cas9 Nickase** and **EnGen Spy dCas9 (SNAP-tag®)** with nuclear localization sequences for genome editing studies. Cas nucleases have been adapted for use in genome engineering, because they can be easily programmed for target specificity by supplying gRNAs of any sequence. In cells and animals, genome targeting is performed by expressing nucleases and gRNA from DNA constructs (plasmid or virus), supplying RNA encoding Cas nuclease and gRNA, or by introducing RNA-programmed Cas nuclease directly.

The latest addition to our EnGen toolbox is **EnGen Lba Cas12a (Cpf1)**, which has a different PAM sequence compared to Cas9. It leaves 5' overhangs after digestion presenting an opportunity for novel CRISPR-based cloning methodologies. It's optimal activity starts already at 16°C making it an extremely efficient choice for researchers working with ectothermic organisms such as zebrafish and xenopus.


EnGen Cas9 Nuclease, *S. pyogenes* (Spy)

Standard genome editing




EnGen Lba Cas12a ("Cpf1")

AT-rich PAM, expanded temperature range




EnGen Spy Cas9 Nickase

increased specificity homology directed repair, dual guide sequence



EnGen Spy dCas9 (SNAP-tag)

in vivo labeling and target enrichment



Available products:

- **EnGen Cas9 NLS, *S. pyogenes***, is an RNA-guided endonuclease that catalyzes site-specific cleavage of double-stranded DNA. The location of the break is within the target sequence 3 bases from the NGG PAM (Protospacer Adjacent Motif).
- **EnGen Spy Cas9 Nickase** is a variant of Cas9 nuclease differing by a point mutation (D10A) in the RuvC nuclease domain, which enables it to nick, but not cleave, DNA. Double-stranded breaks can be generated with reduced off-target cleavage by targeting two sites with EnGen Cas9 Nickase, NLS in close proximity.
- **EnGen Spy dCas9 (SNAP-tag)** is an inactive mutant of Cas9 nuclease that retains programmable DNA binding activity. The N-terminal SNAP-tag allows for covalent attachment of fluorophores, biotin, and a number of other conjugates useful for visualization and target enrichment.
- **EnGen Lba Cas12a (Cpf1)** is a programmable DNA endonuclease guided by a single guide RNA (gRNA). Targeting requires a gRNA complementary to the target site as well as a 5' TTTN protospacer adjacent motif (PAM) on the DNA strand opposite the target sequence. Cleavage by EnGen Lba Cas12a (Cpf1) occurs ~18 bases 3' of the PAM and leaves 5' overhanging ends.

ORDERING INFORMATION

| PRODUCT | NEB # | SIZE |
|-------------------------------------|----------|----------------|
| EnGen® Spy Cas9 Nickase | M0650S/T | 70/400 pmol |
| EnGen® Spy dCas9 (SNAP-tag®) | M0652S/T | 70/400 pmol |
| EnGen® Cas9 NLS, <i>S. pyogenes</i> | M0646T/M | 400/2,000 pmol |
| EnGen® Lba Cas12a (Cpf1) | M0653S/T | 70/2,000 pmol |

Did you get the message?



Everything You Need for Your RNA-Related Workflows

In recent years, the discovery of new classes and modifications of RNA has ushered in a renaissance of RNA-focused research. Did you know that NEB offers a broad portfolio of reagents for the purification, quantitation, detection, synthesis and manipulation of RNA? Experience improved performance and increased yields, enabled by our expertise in enzymology.

Utilize NEB's RNA Toolbox for your RNA related experiments and research and benefit from excellent quality and performance.

Moreover, benefit from our Special Prices for many RNA related products featured in this newsletter. Please ask your local distributor for details.



Did you get the message?
Get special prices on selected NEB RNA Products & Quick Load Purple DNA Ladders!
All eligible discounted products are marked with the "Special Price" icon within this newsletter (pages 10-16). Campaign ends June 30th 2018. Please ask your local distributor for details.

FEATURE ARTICLE – PDF Download

Minding your caps and tails – considerations for functional mRNA synthesis

Applications of synthetic mRNA have grown and become considerably diversified in recent years. Examples include the generation of pluripotent stem cells (1-3), vaccines and therapeutics (4), and CRISPR/Cas9 genome editing applications (5-7). The basic requirements for a functional mRNA – a 7-methylguanylate cap at the 5' end and a poly(A) tail at the 3' end – must be added in order to obtain



efficient translation by eukaryotic cells. Additional considerations can include the incorporation of modified bases, modified cap structures and polyadenylation strategies. Strategies for *in vitro* synthesis of mRNA may also vary according to the desired scale of synthesis. This article discusses options for selection of reagents and the extent to which they influence synthesized mRNA functionality...



Read the full feature article as a PDF on:
www.neb-online.eu/mRNA-Synthesis

Do you want to be on the safe side? Use NEB's murine RNase Inhibitor!

- Specifically inhibits RNases A, B and C.
- Improved resistance to oxidation compared to human/porcine RNase inhibitors
- Ideal for reactions where low DTT concentrations are required (e.g., Real-time PCR)
- Isolated from a recombinant source
- Tested for the absence of DNases and RNases
- No inhibition of polymerase activity used with *Taq* DNA Polymerase, AMV or M-MuLV Reverse Transcriptases, or Phage RNA Polymerases (SP6, T7, or T3)



ORDERING INFORMATION

| PRODUCT | NEB # | SIZE |
|-------------------------|--------|--------------|
| RNase Inhibitor, murine | M0314S | 3,000 units |
| RNase Inhibitor, murine | M0314L | 15,000 units |



RNA Special Price Campaign

High Yield RNA Synthesis in all Flavors

In vitro RNA synthesis requires a DNA template, RNA polymerase, NTPs and other factors. High-yield robust reactions require optimization of each reaction component. NEB offers five *in vitro* RNA synthesis kits, all of which have been optimized to generate reproducible yields of quality RNA.


HiScribe™ T7 High Yield RNA Synthesis Kit & HiScribe™ T7 Quick High Yield RNA Synthesis Kit

The **HiScribe T7 High Yield RNA Synthesis Kit (NEB #E2040)** delivers robust RNA synthesis for a wide range of template sizes. Flexible protocols ensure that performance is maintained even under demanding conditions, such as extended reaction time using very low amounts of template. Protocols are included for partial or complete incorporation of modified or labeled nucleotides in the transcript body, and cap analogs at the RNA 5' end. The **HiScribe T7 Quick High Yield RNA Synthesis Kit (NEB #E2050)** utilizes a master mix format, allowing for faster reaction setup. DNase I and lithium chloride are included for DNA template removal and quick RNA purification.



HiScribe™ T7 ARCA mRNA Synthesis Kit

Most eukaryotic mRNAs require a 7-methyl guanosine (m7G) cap structure at the 5' end and a poly(A) tail at the 3' end for efficient translation to occur. The **HiScribe™ T7 ARCA mRNA Synthesis Kit (NEB #E2060S)** is designed to synthesize capped and tailed mRNAs for a variety of applications. Capped mRNAs are synthesized by co-transcriptional incorporation of Anti-Reverse Cap Analog, ARCA, using T7 RNA Polymerase. A poly(A) tail is then added by *E. coli* Poly(A) Polymerase. **This kit is also available without *E. coli* Poly(A) Polymerase (NEB #E2065S)** for use with DNA templates encoding a poly(A) stretch or not requiring a poly(A) tail. Both kits include DNase I and LiCl for DNA template removal and quick mRNA purification.



Overview of mRNA synthesis work flow with the HiScribe T7 ARCA mRNA Kit (with tailing)

ORDERING INFORMATION


| PRODUCT | NEB # | SIZE |
|------------------------------------------------|--------|---------|
| HiScribe T7 High Yield RNA Synthesis Kit | E2040S | 50 rxns |
| HiScribe T7 Quick High Yield RNA Synthesis Kit | E2050S | 50 rxns |
| HiScribe T7 ARCA mRNA Kit (with Tailing) | E2060S | 20 rxns |
| HiScribe T7 ARCA mRNA Kit | E2065S | 20 rxns |
| HiScribe SP6 RNA Synthesis Kit | E2070S | 50 rxns |



Advantages


- Streamlined format & Quick Workflows
- Flexibility – enables incorporation of cap analogs, radiolabeled and modified nucleotides
- High Yield – up to 180 µg of RNA from a standard 20 µl reaction (e.g. #E2050)
- High Quality Transcripts – optimized formulation for increased transcript integrity
- Get the best translation efficiencies with correctly oriented ARCA caps (#E2060 & #E2065)

Robust RNA Synthesis from a Variety of Template Sizes using the HiScribe T7 High Yield RNA Synthesis Kit



Time course of standard RNA synthesis from three DNA templates of different sizes using HiScribe T7 High Yield RNA Synthesis Kit assembled according to the protocol and incubated at 37°C for the indicated time.

Both cap and tail are required for mRNA function in cell culture.



Purified *Cypridina luciferase* RNA produced as indicated was co-transfected into U2OS cells with purified *Gaussia luciferase* mRNA. mRNAs produced using the HiScribe T7 ARCA mRNA Kit (with Tailing) are 5'-capped and have 3' poly(A) tails. After 16 hours incubation at 37°C, cell culture supernatants from each well were assayed for CLuc and GLuc activity.

Did you get the message?



Monarch Total RNA Miniprep Kit

The Monarch Total RNA Miniprep Kit is a comprehensive solution for sample preservation, cell lysis, gDNA removal, and purification of total RNA from a wide variety of biological samples, including cultured cells, blood, and mammalian tissues. Additionally, tough-to-lyse samples, such as bacteria, yeast, and plant, can be processed with additional steps that enhance lysis. Cleanup of enzymatic reactions or purification of RNA from TRIzol®-extracted samples is also possible using this kit. Purified RNA has high quality metrics, including $A_{260/280}$ and $A_{260/230}$ ratios ≥ 1.8 , high RIN scores, and minimal residual gDNA. Captured RNA ranges in size from full-length rRNAs down to intact miRNAs. Additionally, differential binding conditions allow selective capture or exclusion of the sub-200 nucleotide RNA pool that includes miRNA, 5S rRNA, and tRNA. Purified RNA is suitable for downstream applications, such as RT-qPCR, cDNA synthesis, RNA-seq, Northern blot analysis, etc.


Advantages

- Use with a wide variety of sample types
- Purify RNA of all sizes, including miRNA & small RNAs >20 nucleotides
- Includes DNase I, gDNA removal columns, Proteinase K, and a stabilization reagent
- Protocols available for RNA fractionation and RNA cleanup
- Save money with value pricing for an all-in-one kit

Specifications

- Binding Capacity: 100 μ g RNA
- RNA Size: > 20 nt
- Purity: $A_{260/280}$ and $A_{260/230}$ usually ≥ 1.8
- Input Amount: up to 10^7 or 50 mg tissue*
- Elution Volume: 50 – 100 μ l
- Yield: varies depending on sample type
- Compatible downstream applications: RNA Library prep for NGS, RT-PCR, RT-qPCR, Northern blots
- *View manual for other sample types

Monarch-purified RNA is high-quality and compatible with a wide variety of downstream applications



Total RNA from a broad array of sample types was purified using the Monarch Total RNA Miniprep Kit (NEB #T2010). Aliquots were run on an Agilent® Bioanalyzer® 2100 using the Nano 6000 RNA chip (*S. cerevisiae* RNA was run using a plant Nano assay). RIN values and O.D. ratios confirm the overall integrity and purity of the RNA. To demonstrate compatibility with downstream applications, samples were subsequently used for RT-PCR (+/- RT) (A) for detection of 4 different RNA species using Protoscript® II Reverse Transcriptase (NEB #M0368)/LongAmp® Taq DNA Polymerase (NEB #M0323), NGS library prep (B) using NEBNext® Ultra™ II RNA Library Prep Kit (NEB #E7760) and RT-qPCR (C) using Luna® One-Step RT-qPCR Reagents (NEB #E3005).

| KIT Component | Monarch RNA Purification Kit | Qiagen RNeasy Kits | | |
|------------------------|------------------------------|--------------------|--------------|-----------|
| | | Mini | Protect Mini | Plus Mini |
| gDNA Removal Columns | ✓ | ✗* | ✗* | ✓ |
| DNase I | ✓ | ✗ | ✗ | ✗ |
| Proteinase K | ✓ | ✗ | ✗ | ✗ |
| RNA Protection Reagent | ✓ | ✗ | ✓ | ✗ |

✓ = Included *Not included and not sold separately


Time for change
– try Monarch for free.



For a limited time, try one of our Monarch Nucleic Acid Purification Kits by visiting NEBMonarch.eu


RNA Special Price Campaign

Monarch-purified RNA can be used to prepare high quality RNA-seq libraries for gene expression analysis



Transcript levels in Universal Human Reference RNA (UHRR, Agilent) are compared before and after re-purification using either Qiagen RNeasy® or the Monarch Total RNA Miniprep Kit. Strong correlation with untreated UHRR is observed for both methods (Pearson $R > 0.99$ for both samples). All samples display consistent end-to-end coverage of transcripts indicating an absence of detectable degradation during purification. Poly-A selected RNA was selected from 100 ng of untreated, Qiagen and Monarch samples using the NEBNext Poly(A) mRNA Magnetic Isolation Module (NEB #E7490). RNA-seq libraries were then prepared using the NEBNext Ultra II Directional RNA Library Prep Kit for Illumina® (NEB #E7760) before sequencing on a Miseq® instrument (2 x150). 1.6M reads were randomly sampled from each library and adapter trimmed (Seqprep v1.1). Levels of all Gencode v26 transcripts were assessed using Salmon (0.4) and plotted above (panel A). Average 5'-3' Coverage of Gencode v26 transcripts (assessed by Picard's CollectRnaSeqMetrics 1.56 after mapping to the GRCh38 reference genome with Hisat v2.0.3 and marking duplicates with Picard's MarkDuplicates 1.56) is shown below (panel B).

The Monarch Total RNA Miniprep Kit successfully purifies small RNAs below 200 nucleotides, enabling a more accurate representation of the total RNA pool



RNA preps were performed on HEK293 cells, mouse heart, rat brain, or rat spleen using the Monarch Total RNA Miniprep Kit (N) (NEB #T2010) and the RNeasy® Mini Kit from Qiagen (Q). Equivalent amounts were resolved on a Bioanalyzer 2100 using the Small RNA chip. Monarch-purified RNA contains significantly more RNA in the sub-200 nucleotide pool.



TIPS FOR SUCCESSFUL RNA EXTRACTIONS

- **Inactivate RNases after harvesting your sample:** Nucleases in your sample will lead to degradation, so inactivating them is essential. Process samples quickly, or use preservation reagents, and always ensure nuclease-free working environments.
- **Do not exceed recommended input amounts:** Buffer volumes are optimized for the recommended input amounts. Exceeding these can result in inefficient lysis and can also clog the column.
- **Ensure samples are properly homogenized/ disrupted:** Samples should be disrupted and homogenized completely to release all RNA.
- **For sensitive applications, ensure proper gDNA removal:** gDNA removal is removed by the gDNA removal column and subsequent on-column DNase I treatment. Off-column DNase I treatment can also be employed.



Request a free sample*!
Visit nebmonarch.eu
or contact your local distributor!



ORDERING INFORMATION

| RNA PURIFICATION KIT | NEB # | SIZE |
|--------------------------------|--------|----------|
| Monarch Total RNA Miniprep Kit | T2010S | 50 preps |

COLUMNS AVAILABLE SEPARATELY

| | | |
|----------------------------------|--------|-----------------------|
| Monarch RNA Purification Columns | T2007L | 100 columns and tubes |
| Monarch gDNA Removal Columns | T2017L | 100 columns and tubes |
| Monarch Collection Tubes II | T2018L | 100 tubes |

BUFFERS & REAGENTS AVAILABLE SEPARATELY

| | | |
|-------------------------------------------------------------------------------------------|----------|--------------|
| Monarch DNA/RNA Protection Reagent | T2011L | 56 ml |
| Monarch RNA Lysis Buffer | T2012L | 100 ml |
| Monarch Total RNA Miniprep Enzyme Pack (contains DNase I, Prot K, and associated buffers) | T2019L | 1 pack |
| Monarch RNA Priming Buffer | T2013L | 56 ml |
| Monarch RNA Wash Buffer | T2014L | 50 ml |
| Nuclease-free Water | B1500S/L | 25 ml/100 ml |

NEB also offers an extensive range of Monarch DNA Purification Kits and products. Visit www.nebmonarch.eu to learn more!


Did you get the message?



RNA Quantitation by RT-qPCR: Novel enzyme technology for superior RT-qPCR Performance

New England Biolabs offers a bright new choice for your qPCR and RT-qPCR. Luna products have been optimized for robust performance on diverse sample sources and target types. Available for dye-based or probe-based detection, Luna products can be used across a wide variety of instrument platforms.

For superior RT-qPCR Performance, experience a novel, more thermostable "Designer" Reverse Transcriptase ("WarmStart RT") as an integral part of the Luna Universal One-Step RT-qPCR Kits. This new technology highly improves the performance and is paired with Hot Start *Taq* which increases reaction specificity and robustness.



A blue visible dye assists in tracking the reagents when pipetting into clear, multi-welled PCR plates.

EXPERIENCE BEST-IN-CLASS PERFORMANCE

- All Luna products have undergone rigorous testing to optimize specificity, sensitivity, accuracy and reproducibility.
- Products perform consistently across a wide variety of sample sources.
- A comprehensive evaluation of commercially-available qPCR and RT-qPCR reagents demonstrates superior performance of Luna products.

OPTIMIZE YOUR RT-qPCR WITH LUNA WARMSTART® REVERSE TRANSCRIPTASE

- Novel, thermostable reverse transcriptase (RT) improves performance.
- WarmStart RT paired with Hot Start *Taq* increases reaction specificity and robustness.
- Convenient and fast One-Step RT-qPCR set-up reducing the risk of pipetting errors


MAKE A SIMPLER CHOICE

- One product per application simplifies selection.
- Convenient master mix formats and user-friendly protocols simplify reaction setup.
- Non-interfering, visible tracking dye helps to eliminate pipetting errors.

Find the right Luna product for your application:

| 1 Select your target | 2 Select your detection method | |
|----------------------|--------------------------------------------------|--------------------------------------------------------|
| | Dye-based | Probe-based |
| Genomic DNA or cDNA | Luna Universal qPCR Master Mix (NEB #M3003) | Luna Universal Probe qPCR Master Mix (NEB #M3004) |
| RNA | Luna Universal One-Step RT-qPCR Kit (NEB #E3005) | Luna Universal Probe One-Step RT-qPCR Kit (NEB #E3006) |

Luna Universal One-Step RT-qPCR Kit offers exceptional sensitivity, reproducibility & RT-qPCR performance



RT-qPCR targeting human GAPDH was performed using the Luna Universal One-Step RT-qPCR Kit (Input: 1 µg – 0.1 pg Jurkat total RNA); NTC = non-template control



Request a free sample*!
Visit LUNAqPCR.eu
or contact your local distributor



ORDERING INFORMATION

| PRODUCT | NEB # | SIZE |
|-------------------------------------------|--------------|--------------------------|
| Luna Universal qPCR Master Mix | M3003S/L/X/E | 200/500/1,000/2,500 rxns |
| Luna Universal Probe qPCR Master Mix | M3004S/L/X/E | 200/500/1,000/2,500 rxns |
| Luna Universal One-Step RT-qPCR Kit | E3005S/L/X/E | 200/500/1,000/2,500 rxns |
| Luna Universal Probe One-Step RT-qPCR Kit | E3006S/L/X/E | 200/500/1,000/2,500 rxns |


*Limit one per customer, while supplies last.

RNA Special Price Campaign

Affordable Performance cDNA Synthesis: ProtoScript® II RTase & cDNA Kits


ProtoScript II Reverse Transcriptase is a recombinant M-MuLV reverse transcriptase with reduced RNase H activity and increased thermostability. It can be used to synthesize first strand cDNA at higher temperatures than the wild-type M-MuLV. The enzyme is active up to 48°C, providing higher specificity, higher yield of cDNA and more full-length cDNA product, up to 12 kb in length.

Take advantage of the low cost per unit with ProtoScript II




*based on published list prices as of 01/2018; prices may differ in various European countries!

cDNA synthesis for RT-PCR up to 10 kb using the ProtoScript II First Strand cDNA Synthesis Kit



Jurkat total RNA (250 ng) was converted to first strand cDNA using the ProtoScript II First Strand cDNA Synthesis Kit. Using 10% of the first strand cDNA reactions, amplicons representing four different mRNAs were amplified using LongAmp® Taq 2X Master Mix (NEB #M0287). Lane 1, (mRNA) SDHA, 1.9 kb; lane 2, (mRNA) HERC1, 5.5 kb; lane 3, (mRNA) XRN1, 7.3 kb; lane 4, (mRNA) FBN1, 9.2 kb; Ladder (L) is a 2-log DNA ladder (NEB #N3200).



Jurkat total RNA (1 µg) was used in a 20 µl first strand cDNA synthesis. Mixtures of all reaction components, except for reverse transcriptase, were held at different temperatures for 3 min. 200 units NEB's ProtoScript® II Reverse Transcriptase (A) or SuperScript® II (B) was added and incubated at the indicated temperature for 50 minutes, followed by heat inactivation for 5 min at 80°C. 1 µl of cDNA was used in a 25 µl PCR using LongAmp® Hot Start Taq 2X Master Mix (NEB #M0533) for 35–40 cycles. Ladder M is the Quick-Load® 2-Log DNA Ladder (NEB #N0469).

Advantages

- High cDNA yield
- Superior performance for longer templates (up to 12 kb)
- Increased thermostability (up to 50°C)
- Value pricing
- Also available as affordable first strand cDNA synthesis kit

TIPS FOR SUCCESSFUL cDNA Synthesis

- You can increase the yield of long cDNA products by doubling the amount of enzyme and dNTP.
- Intact RNA of high purity is essential for generating cDNA for cloning applications.
- Total RNA or mRNA can be used in the reverse transcription reaction. Total RNA is generally sufficient for cDNA synthesis reactions. However, if desired, mRNA can be easily obtained using a polyA Spin mRNA Isolation Kit (NEB #S1560) or Magnetic mRNA Isolation Kit (NEB #S1550).
- The amount of RNA required for cDNA cloning depends on the abundance of the transcript-of-interest. In general 1 ng to 1 µg total RNA or 0.1 ng to 100 ng mRNA are recommended.


ORDERING INFORMATION

| PRODUCT | NEB # | SIZE |
|-------------------------------------------------|------------|---------------------------|
| ProtoScript® II First Strand cDNA Synthesis Kit | E6560S/L | 30/150 rxns |
| ProtoScript® II Reverse Transcriptase | M0368S/L/X | 4,000/10,000/40,000 units |



Do you use DNA ladders in your daily lab work?

We thought so! We're confident that you will love our Quick-Load[®] purple DNA ladders, with their sharp bands, ready-to-use format, and competitive pricing. In addition, our purple ladders cast no UV shadow, so you'll never miss a band.



1 kb DNA Ladder

0.8% TAE agarose gel.
Mass values are for
0.5 µg/lane.

NEB #N0552

100 bp DNA Ladder

1.3% TAE agarose gel.
Mass values are for
0.5 µg/lane.

NEB #N0551

2-Log DNA Ladder

1.0% TBE agarose gel.
Mass values are for
1 µg/lane.

NEB #N0550

50 bp DNA Ladder

3% TBE agarose gel.
Mass values are for
1 µg/lane.

NEB #N0556

Low Molecular Weight DNA Ladder

3% TBE agarose gel.
Mass values are for 0.5
µg/lane.

NEB #N0557



UV shadow comparison

The innovative Gel Loading Dye, Purple (6X) (#B7024S) (Lane 1) included in the Quick-Load Purple 100 bp DNA Ladder does not cast a UV shadow over the underlying bands, unlike conventional Bromophenol-Blue containing gel loading dyes (Lane 2).

← UV shadow

ORDERING INFORMATION

| PRODUCT | NEB # | SIZE |
|----------------------------------------------------------------|----------|-------------------|
| Quick-Load [®] Purple 2-Log DNA Ladder | N0550S/L | 250/750 gel lanes |
| Quick-Load [®] Purple 100 bp DNA Ladder | N0551S/L | 125/375 gel lanes |
| Quick-Load [®] Purple 1 kb DNA Ladder | N0552S/L | 125/375 gel lanes |
| Quick-Load [®] Purple 50 bp DNA Ladder | N0556S | 250 gel lanes |
| Quick-Load [®] Purple Low Molecular Weight DNA Ladder | N0557S | 125 gel lanes |

SPECIAL PRICE*
until 30.06.2018

Your local NEB distributor:

BeNeLux:
BIOKÉ
Tel: (+31) 71 720 0220
(BE: 0800 - 71640)
Fax: (+31) 71 891 0019
info@bioke.com
www.bioke.com

BULGARIA
ELTA 90M Ltd.
Tel: (+359) 2 983 96 49
Fax: (+359) 2 9832211
office@elta90.eu
www.elta90.com

CROATIA:
DIAGNOSTICA SKALPELI D.O.O.
Free call: 800-315911
Tel: (+385) 1 377 8484
Fax: (+385) 1 377 8585
skalpeli@skalpeli.hr
www.skalpeli.hr

CZECH REPUBLIC:
BIOTECH A.S.
Tel: (Toll Free) 0800 124683
Fax: 02 72701742
info@ibiotech.cz
www.ibiotech.cz

DENMARK:
BioNordika Denmark A/S
DENMARK
Tel: (39) 56 20 00
Fax: (39) 56 19 42
info@bionordika.dk
www.bionordika.dk

FINLAND:
BioNordika Oy
Tel: (+358) 207 410 270
info@bionordika.fi
www.bionordika.fi

GREECE:
BIOLINE SCIENTIFIC
Tel: (0030) 210 522 6547
Fax: (0030) 210 524 4744
info@bioline.gr
www.bioline.gr

HUNGARY:
KVALITEX SCIENTIFIC
Phone: (1) 340 4700
Fax: (1) 339 8274
info@kvalitex.hu
www.kvalitex.hu

ITALY:
EUROCLONE S.P.A.
Free call: 800-315911
Tel: +39 02 38195.1
Fax: +39 02 3810145
info@euroclone.it
www.euroclonegroup.it

ISRAEL:
ORNAT LTD.
Tel: (+972) 8 947 7077
Fax: (+972) 8 936 3034
ornatbio@ornat.co.il
www.ornat.co.il

NORWAY:
BioNordika AS
Tel +47 23 03 58 00
Fax +47 23 03 58 01
info@bionordika.no
www.bionordika.no

POLAND:
Lab-JOT Ltd. Sp.z o.o. Sp.k.
Stacjonarny: +48 22 335 988 4
Komórkowy: +48 606 338 879
Fax: +48 22 335 981 9
labjot@labjot.com
www.labjot.com

PORTUGAL:
Werfen
Tel. +351 214 247 312
Fax +351 21 425 52 80
dac-pt@werfen.com
http://pt.werfen.com

ROMANIA:
ELTA 90 MEDICAL RESEARCH SRL
Tel: +40 21 2322694
Fax: +40 21 2322696
office@elta90mr.ro
www.elta90mr.ro

RUSSIA:
SkyGen LLC
Tel: +7 495 215 02 22
info@skygen.com
www.skygen.com

SERBIA:
ALFA GENETICS D.O.O.
Tel: 00381 11 2084252
Fax: 00381 11 2084253
office@alfagenetics.rs
www.alfagenetics.rs

SPAIN:
Werfen
Tel: +34 902 20 30 90
Fax: +34 902 22 33 66
customerservice-es@werfen.com
http://es.werfen.com

SLOVENIA:
Mikro+Polo d.o.o.
Tel: +386 2 614 33 00
Fax: +386 2 614 33 20
info@mikro-polo.si
www.mikro-polo.com

SOUTH AFRICA:
INQABA BIOTEC
Tel: +27 123 43 58 29
Fax: +27 86 677 8409
info@inqababiotec.co.za
www.inqababiotec.co.za

SWEDEN, ESTONIA, LATVIA,
LITHUANIA:
BIONORDIKA SWEDEN AB
Tel: +46 8 30 60 10
Fax: +46 8 30 60 15
info@bionordika.se
www.bionordika.se

SWITZERLAND:
BIOCONCEPT
Tel: +41 (0)61 486 80 80
Fax: +41 (0)61 486 80 00
info@bioconcept.ch
www.bioconcept.ch

TURKEY:
SACEM HAYAT TEKNOLOJILERI
Gebze Plastikçiler Organize Sanayi
Bölgesi (GEPOSB)
Tel: 0262 751 02 74
Fax: 0262 751 02 75
sacem@sacem.com.tr
www.sacem.com.tr



international.neb.com

