

NEB expressions

a scientific update



in this issue

- 2** Substrate specificity and mismatch discrimination in DNA ligases
- 5** HiFi *Taq* DNA Ligase
- 6** New products for isothermal amplification
- 8** Tools to support genome editing workflows
- 10** Lighting the way: New LUNA universal qPCR & RT-qPCR Reagents
- 12** Request your free LUNA qPCR or RT-qPCR Sample

Substrate specificity and mismatch discrimination in DNA ligases

DNA ligases vary in their ability to join fragments, add adaptors, repair nicks and breaks, link vectors and inserts, and to circularize dsDNA. Ligases also vary in their type of activity. The specificity and accuracy of the ligation depends upon ligase selection and careful optimization of reaction conditions. With the right ligase, conditions and probes, even single-base variations in sequence can be reliably detected.

DNA ligases are enzymes that seal breaks in DNA by joining 5'-phosphorylated DNA termini to 3'-OH DNA termini (1-4). *In vivo*, ligases are important for the repair of nicks, breaks in one strand of a dsDNA molecule, DNA formed during replication (i.e., Okazaki fragments), as well as both nick and double-strand break joining during repair events (5). *In vitro*, ligases (notably T4 DNA Ligase) are critical reagents for many molecular biology protocols, including vector-insert joining for recombinant plasmid construction, adaptor ligation for next-generation sequencing (NGS) library construction, and circularization of dsDNA (6). T4 DNA Ligase (NEB #M0202) is incredibly efficient at sealing nicks, as well as joining or circularizing DNA fragments with blunt or cohesive (short complementary) ends. This activity can be further improved with the addition of macromolecular enhancers, such as polyethylene glycol (PEG), as seen in NEB's Quick Ligation™ Kit (NEB #M2200) (7,8).

Less commonly utilized *in vitro*, *Taq* DNA Ligase (NEB #M0208) will ligate only nicks (9-12). *Taq* Ligase is a NAD⁺-dependent DNA ligase from a thermostable bacterium that can survive high temperatures (up to 95°C) and is active over a range of elevated temperatures (37–75°C). However, it only has significant activity on nicked DNA, and negligible activity on short cohesive and blunt substrates in end-joining reactions. Given these limitations, and the fact that T4 DNA Ligase can ligate everything *Taq* ligates, and many more structures, why not use T4 DNA Ligase for all applications?

T4 DNA Ligase can ligate a wide variety of DNA structures, including modified bases and the ends of double stranded fragments. It will also efficiently ligate many undesirable structures, including substrates containing gaps of one or more nucleotides and nicked substrates that contain DNA base pair mismatches (12-15). In most cases, this unwanted activity isn't a problem,

for example, when joining 1 or 2 fragments into a plasmid, or pushing an adaptor ligation reaction as far towards completion as possible to prepare high yields of DNA NGS libraries.

For some applications, however, there cannot be any end-joining activity at all, and for others, there is a need for the exclusive ligation of fully base-paired nicks with no gaps. For example, DNA assembly methods, such as Gibson Assembly® (NEB #E5510) and NEBuilder® HiFi DNA Assembly (NEB #E2621), require nick-selective ligases; this method utilizes long overlaps that are dynamically generated by exonucleases and gaps are filled by a DNA polymerase (16,17). Final joining is accomplished by a nick-selective ligase, such as *Taq* DNA Ligase, which only reacts with substrates containing no gaps, and will not join any fragments end-to-end without the exo/polymerase generation of annealed complementary regions. The use of a nick-selective ligase ensures that fragments are not joined out of order, and no deletions result from ligation across nucleotide gaps in annealed structures. (For more information, see www.neb.com/DNAassembly).


LIGASE SPECIFICITY

DNA ligases generally prefer fully Watson-Crick base-paired dsDNA substrates to those containing one or more mismatches. However, ligases can ligate some mismatches to a significant degree, and very active ligases, such as T4 DNA Ligase, can ligate nicks containing one or more mismatches near the ligation junction with high efficiency (15,18). Ligases are thought to interrogate dsDNA for proper base pairing through minor groove contacts, and thus do not read specific base sequences, but are sensitive to distortions in helix shape (19). Large purine:purine mismatches and most smaller pyrimidine:pyrimidine mismatches are typically worse ligation substrates than pyrimidine:purine mismatches. Helix stability also plays some role, and mismatches with more hydrogen bonds are more readily ligated than those with few. For many ligases, G:T mismatches, with two hydrogen bonds and a base-pair size nearly indistinguishable from a Watson-Crick base pair, are joined with nearly the same efficiency as a correct base pair. Additionally, DNA ligases have been generally found to have a higher discrimination at the upstream side of the ligation junction (the base pair providing the 3'-OH terminus to the ligation) than on the downstream side (the base pair providing the 5'-phosphate to the ligation). The structural/mechanistic reason for this differential is not known for certain, but may have to do with the slight melting of the 5'-terminus during the reaction. This "peeling back" of the 5'-phosphorylated base can be observed in the crystal structures of several DNA ligases bound to substrate (20,21).

Thermostable DNA ligases, including *Taq* DNA ligase, are naturally able to discriminate against


 **FIGURE 1: Ligase Detection Reaction (LDR)**

Two ligation probes are designed such that they are complementary to a target region of interest, and anneal with no gaps. Typically, if a SNP is to be resolved, the nucleotide of interest is situated at the junction of the two probes. The probes are combined with the DNA to be examined (typically genomic or PCR amplified region) and the thermostable DNA ligase. The DNA is melted, then cooled to a ligation temperature that allows the probes to anneal to the target. If the probes anneal to form a nicked structure with no gaps or mismatches, efficient ligation will proceed. Cycling, melting and annealing/ligation allows successive rounds of probes to anneal and ligate, resulting in linear amplification of the ligation product if the sequence of interest is present in the target DNA.



**FIGURE 2:
Ligase Chain Reaction (LCR)**

In this method, two pairs of probes are designed, one pair complementary to the top strand of the target, one pair to the bottom. Upon melting and annealing, both probe pairs can anneal to the target, and ligate efficiently if they form a nicked sequence without gaps or mismatches. On successive rounds of melting and re-annealing, unligated probe can now anneal to both the original target DNA and to the probes ligated in previous rounds. As each ligation product becomes a template for the complementary probe pair, LCR enables exponential amplification of the ligated product.



ligating substrates containing base pair mismatches (i.e., are "higher fidelity") than T4 DNA Ligase (18,22,23). Despite this higher fidelity, *Taq* DNA Ligase can still detectably ligate many T:G, T:T, and A:C mismatches. Thermostable DNA ligases are active at elevated temperatures, allowing further discrimination by incubating the ligation at a temperature near the melting temperature (T_m) of the DNA strands. This selectively reduces the concentration of annealed mismatched substrates (expected to have a slightly lower T_m around the mismatch) over annealed fully base-paired substrates. Thus, high-fidelity ligation can be achieved through a combination of the intrinsic selectivity of the ligase active site and careful balance of reaction conditions to reduce the incidence of annealed mismatched dsDNA.


APPLICATIONS REQUIRING HIGH-FIDELITY LIGATION

Numerous applications have been developed that take advantage of the high fidelity of *Taq* and other thermostable DNA ligases to detect specific nucleotide sequences with high specificity and quantitative accuracy, including profiling single nucleotide polymorphisms (SNPs) (9,24,25). In the Ligase Detection Reaction (LDR), a set of probes complementary to the sequence of interest are annealed to target DNA (genomic DNA, or a PCR amplified fragment) and treated with a high-fidelity thermostable DNA ligase (Figure 1, page 2). If the target sequence is present, the probes will ligate; cycling through rounds of melting and annealing can allow linear

amplification of the probe ligation product. With the right ligase, conditions and suitable probes, single-base differences can be reliably detected. The original paper detected the ligation product through visualization in a gel, but detection through fluorophore-quencher pairs or qPCR-based methods can greatly increase the sensitivity of detection (26-32). LDR has also been

**FIGURE 3:
Ligation-Rolling Circle Amplification/Padlock Probes**

In this method, a single probe is designed such that the ends of the probes are complementary to the target sequence. When annealed to the desired target, the ends form a nicked structure that can be efficiently ligated if there are no gaps or mismatched base pairs. Exonuclease treatment destroys the uncircularized DNA, and the remaining circular structures can be detected through rolling circle amplification (RCA), or linearized and amplified with PCR.



extended to multiplexed probe sets that allow the simultaneous interrogation of multiple potential SNP sites (27).


The closely related Ligase Chain Reaction (LCR) takes the LDR method and makes it amplifiable in an exponential fashion (9,24). In LCR, four probes are used, one pair complementary to one target strand, and a second pair of probes complementary to the other strand (Figure 2). Since the probe pairs are complementary to each other, the probe ligated in one cycle becomes a template for ligation of additional probe in subsequent cycles. This methodology allows for detection of SNPs with greater sensitivity than the original LDR method, but requires extraordinary discrimination against mismatch ligation for both probe sets, as even trace ligation on a mismatched template will result in template for further probe amplification (and thus, an erroneous positive signal). The complementarity of LCR probe pairs also means that probes can anneal to each other, forming blunt-end or single-base overhang substrates, depending on probe design strategy. While most high-fidelity ligases have far lower activity on double-stranded fragment end joining, even trace blunt-end activity will generate template for further rounds of high-efficiency nick ligation. Thus, LCR can suffer from high non-templated background as well, and requires careful probe design. The modification gap-LCR method attempts to address these background and discrimination issues by utilizing probes which anneal with a single-nucleotide gap that must be filled by a polymerase in order to generate a substrate suitable for ligation (33). This modification leverages the discrimination against gap ligation of thermostable high-fidelity ligases, but requires a thermostable polymerase and dNTPs as well.



FIGURE 4:

Comparison of fidelity of *Taq* DNA Ligase (NEB #M0208), Ampligase (Epicentre), and HiFi *Taq* DNA Ligase (NEB #M0647)

Fidelity measurements were performed using 1 μ L of ligase in a 50 μ L reaction mixture in the supplied buffers at 1x concentration. Reactions were incubated 30 min at 55°C, using multiplexed substrate pools as outlined in our previous publication (Lohman, G.J. et al. 2016). The rows represent a single template sequence, while columns indicate a particular ligation product resulting from a specific pair of probes ligating with the indicated bases at the ligation junction. A dot indicates detection of a product (see legend). The diagonal from the top left to the bottom right represents Watson-Crick ligation products; all other spaces indicate mismatch ligation products. While *Taq* ligase and Ampligase perform similarly under these conditions, with a range of mismatch products detectable, HiFi *Taq* Ligase shows dramatically fewer mismatch products while maintaining high yields (image adapted from Lohman, G.J. et al. 2016).



Additional detection-by-ligation technologies have been devised to take advantage of high-fidelity ligation events by generating circular templates that can be detected in a secondary reaction (34-38). In the “padlock” probe design, a single-stranded probe is devised where the 5’ and 3’ ends are both complementary to a target sequence (Figure 3, page 3). Much like LDR, the ends of this single probe form a nick structure with no gaps when annealed, and can only be ligated when fully base-paired to a complementary target. All single-stranded DNA can be destroyed by exonuclease treatment, and the circularized probes can be detected by methods such as rolling circle amplification, or linearized and detected by PCR. Variation on the padlock probe design include “molecular inversion probes,” which, similar to gap-LCR, have a single-stranded probe in which both ends are annealed with a gap of one or more nucleotides that must be filled by a polymerase before ligation can occur.

OTHER FACTORS TO CONSIDER

These and other detection-by-ligation methodologies depend on the ability of the ligase to discriminate against substrates containing one or more mismatches, yet retain high activity on even low concentrations of the fully base-paired probe-target structure. While the choice of ligase is very important, careful probe design, selection of reaction temperature, and even ligation buffer conditions can all contribute to the fidelity of the ligation reaction, and thus the accuracy and sensitivity of the detection-by-ligation. For example, probes should take advantage of the naturally higher discrimination of ligases on the

upstream side of the ligation junction (the base pair providing the 3’-hydroxyl). Probes that place the base of interest on the downstream side will provide significantly poorer discrimination, and probes will ligate on templates containing other bases at the targeted position. Furthermore, it is important to know what base pair mismatches are more easily ligated by a given ligase. For example, if you are targeting a position that can be an A:T or a T:A base pair, it would be better to use a probe with an A at the 3’ end (targeting the strand with a T at the SNP position) than to use a probe with a T; when annealing to the wrong SNP, the first case would result in a difficult-to-ligate A:A mismatch, while the second would result in a T:T mismatch that can be ligated by *Taq* DNA Ligase with relatively high efficiency.

Incubation temperature is also a key consideration, and typically must be optimized for each application. If the ligation temperature is too far below the T_m of the probes, even mismatched probes will be annealed, increasing the chances of ligation occurring. If the ligation temperature is too far above the T_m, fully complementary sequences will not be annealed. High-fidelity ligation reactions should typically be run 1-2°C below the T_m of the probes to give the highest possible accuracy by minimizing the concentration of annealed mismatched probes. Consequently, it is important to match the T_m of both the upstream and downstream probe annealing regions, and all probe sets when attempting a multiplexed reaction. If there is a range of T_m values for the probe annealing regions, no single reaction temperature will result in the optimal balance of fidelity versus activity for all probe sets.

Buffer conditions can also affect the fidelity of DNA ligases. In particular, it has been observed for several thermostable DNA ligases, including T4 DNA Ligase and Human DNA Ligase 3, that increasing monovalent cation concentration improves the fidelity of ligation (15,39). This effect is thought to be related to a weakening of the binding of the ligase to its substrate, with a disproportionate suppression binding/ligation of mismatched substrates. Too much salt can erase activity on even fully base-paired substrates, and the best salt balance for each ligase must be empirically determined.

OPTIMIZATION THROUGH HIGH-THROUGHPUT PROFILING

NEB researchers recently published a method for the high-throughput profiling of ligase fidelity, a method that extends earlier studies through a high-sensitivity multiplexed format (18,40). This methodology has been used to rapidly screen buffer conditions for *Taq* DNA Ligase and their effect on fidelity. In this method, substrate pools were prepared consisting of one target (template) strand and four upstream probes and four downstream probes, each differing only in the base at the ligation junction. Thus, all four bases at either side of the ligation junction were represented. Sixteen separate pools were prepared, each with a different template strand covering all 16 possible NN pairs in the template as well. The probes were designed such that each possible pairing resulted in a product of unique length, with products repeatable and quantifiable by capillary electrophoresis (CE). This method allowed screening of all possible base combina-

tions (Watson-Crick and mismatched) around the ligation junction in 16 wells of a 96-well plate, allowing 6 conditions to be screened per plate. The results indicated that the optimal buffer for *Taq* DNA Ligase contains with 100 – 200 mM KCl at pH 8.5.

CONCLUSION

This optimization method has been used internally at NEB to screen additional ligases, conditions and formulations, and has led to the development of the new HiFi *Taq* DNA Ligase (NEB #M0647, see below). Using this method, both the enzyme and the reaction buffer were optimized, resulting in the highest fidelity NAD⁺-dependent DNA ligase commercially available (Figure 4).

It is important to note that thermostable, high-fidelity, nick-selective DNA ligases like *Taq* DNA Ligase, HiFi *Taq* DNA Ligase, and the ATP dependent 9°N™ DNA Ligase (NEB #M0238), are not replacements for T4 DNA Ligase in applications such as routine cloning or DNA library preparation. However, when a method relies on accurate ligation of nicks lacking gaps

or mismatched base pairs, using one of these ligases, combined with careful probe design and reaction condition optimization, will be critical for success.

References

- Shuman, S. and C. D. Lima (2004). *Curr. Opin. Struct. Biol.* 14, 757-764.
- Tomkinson, A. E., et al. (2006). *Chem. Rev.* 106, 687-699.
- Ehlenberger, T. and A. E. Tomkinson (2008). *Annu. Rev. Biochem.* 77, 313-338.
- Shuman, S. (2009). *J. Biol. Chem.* 284, 17365-17369.
- Arakawa, H. and G. Iliakis (2015). *Genes (Basel)* 6, 385-398.
- Lohman, G. J., et al. (2011). *Curr. Protoc. Mol. Biol.* Chapter 3: Unit 3.14.
- Pheiffer, B. H. and S. B. Zimmerman (1983). *Nucleic Acids Res.* 11, 7853-7871.
- Hayashi, K., et al. (1986). *Nucleic Acids Res.* 14, 7617-7631.
- Barany, F. (1991). *Proc. Natl. Acad. Sci. USA* 88, 189-193.
- Barany, F. and D. H. Gelfand (1991). *Gene* 109, 1-11.
- Lauer, G., et al. (1991). *J. Bacteriol.* 173, 5047-5053.
- Luo, J. and F. Barany (1996). *Nucleic Acids Res.* 24, 3079-3085.
- Nilsson, S. V. and G. Magnusson (1982). *Nucleic Acids Res.* 10, 1425-1437.
- Goffin, C., et al. (1987). *Nucleic Acids Res.* 15, 8755-8771.
- Wu, D. Y. and R. B. Wallace (1989). *Gene* 76, 245-254.
- Gibson, D. G., et al. (2009). *Nat Methods* 6, 343-345.
- Gibson, D. G., et al. (2010). *Science* 329, 52-56.
- Lohman, G. J., et al. (2016). *Nucleic Acids Res.* 44, e14.
- Liu, P., et al. (2004). *Nucleic Acids Res.* 32, 4503-4511.
- Pascal, J. M., et al. (2004). *Nature* 432, 473-478.
- Nandakumar, J., et al. (2007). *Mol. Cell* 26, 257-271.
- Luo, J., et al. (1996). *Nucleic Acids Res.* 24, 3071-3078.
- Nishida, H., et al. (2005). *Acta. Crystallogr. Sect. F Struct. Biol. Cryst. Commun.* 61, 1100-1102.
- Barany, F. (1991). *PCR Methods Appl.* 1, 5-16.
- Wiedmann, M., et al. (1994). *PCR Methods and Applications* 3, S51-64.
- Zirvi, M., et al. (1999). *Nucleic Acids Res.* 27, e40.
- Pingle, M. R., et al. (2007). *J. Clin. Microbiol.* 45, 1927-1935.
- Cheng, C., et al. (2013). *Anal. Biochem.* 434, 34-38.
- Hamada, M., et al. (2013). *Electrophoresis* 34, 1415-1422.
- Hommatsu, M., et al. (2013). *Anal. Sci.* 29, 689-695.
- LeClair, N. P., et al. (2013). *J. Clin. Microbiol.* 51, 2564-2570.
- Watanabe, S., et al. (2014). *Anal. Chem.* 86, 900-906.
- Marshall, R. L., et al. (1994). *PCR Methods and Applications* 4, 80-84.
- Nilsson, M., et al. (1994). *Science* 265, 2085-2088.
- Cao, W. (2001). *Clin. Appl. Immun. Rev.* 2, 33-43.
- Qi, X., et al. (2001). *Nucleic Acids Res.* 29, E116.
- Cao, W. (2004). *Trends in Biotechnology* 22, 38-44.
- Cheng, Y., et al. (2013). *Analyst* 138, 2958-2963.
- Bhagwat, A. S., et al. (1999). *Nucleic Acids Res.* 27, 4028-4033.
- Greenough, L., et al. (2016). *Nucleic Acids Res.* 44, e15.

Need help calculating your optimal ligation temperature?

Try our
Thermostable Ligase
Reaction Temp. Calculator
at ligasecalc.neb.com




NEW PRODUCT

HiFi *Taq* DNA Ligase

Optimized HiFi *Taq* DNA Ligase efficiently seals nicks in DNA with unmatched fidelity (See Figure 4, page 4). Ligation of two adjacent oligonucleotides that are hybridized to a complementary target DNA is enhanced by an improved reaction buffer, and mismatch ligation is dramatically reduced. Furthermore, HiFi *Taq* DNA Ligase exhibits increased thermostability and is active at elevated temperatures (37–75°C) for extended cycles. This improved formulation allows for higher resolution discrimination between ligation donors and acceptors at both the 3' and 5' side of the ligation junction, enabling precise detection of SNPs and other allele variants via detection methods such as LCR and LDR.

HiFi *Taq* DNA Ligase exhibits increased thermostability

HiFi *Taq* DNA Ligase and Ampligase® (1 µL enzyme in a 50 µL reaction) were cycled (80°C for 90 seconds/94°C for 10 seconds) up to 100 times in their respective 1X reaction buffer. Ligase activity was assayed using a FAM-labeled nicked dsDNA substrate detected by capillary electrophoresis.



advantages

- High fidelity
- Increased thermostability (as compared to Ampligase®)
- Recombinant source
- Improved discrimination at both the 3' and 5' side of the ligation junction
- Online tool to calculate optimal ligation temperature (ligasecalc.neb.com)

ORDERING INFORMATION:

PRODUCT	NEB #	SIZE
HiFi <i>Taq</i> DNA Ligase	M0647S	50 rxns

Watch the video: www.neb-online.eu/ligase




New Products for Isothermal Amplification

WarmStart® LAMP Kit (DNA & RNA)

Loop Mediated Isothermal Amplification (LAMP) is a commonly-used technique for rapid nucleic acid detection. NEB's WarmStart LAMP products provide a simple, one-step solution for DNA or RNA targets. The master mix supplied with the WarmStart LAMP Kit contains the robust and rapid Bst 2.0 WarmStart DNA Polymerase and WarmStart RTx Reverse Transcriptase, both *in silico*-designed enzymes for improved performance in LAMP reactions. The kit also includes a fluorescent dye to enable real-time fluorescence measurement of LAMP. The WarmStart LAMP Kit is compatible with multiple detection methods, as shown in the figure below.


 NEB's WarmStart LAMP Kit (DNA & RNA) is compatible with multiple detection methods*



looking for best-in-class visual detection?


Try our WarmStart Colorimetric LAMP 2X Master Mix (DNA & RNA)

- NEB's WarmStart Colorimetric LAMP 2X Master Mix (DNA & RNA) offers the same robust performance as the WarmStart LAMP Kit, and contains a colorimetric dye for best in class visual detection of your target.




Yellow color indicates that amplification has occurred.

Turbidity detection



LAMP product is visible by eye due to the precipitation of magnesium pyrophosphate

Endpoint visualization




Visualization using a fluorescent intercalating dye (e.g., ethidium bromide, SYBR® green or calcein/HNBlue)

* The NEB WarmStart LAMP Kit (DNA & RNA) includes separate fluorescent dye for real-time fluorescence measurement. Alternately, detection can be accomplished by turbidity detection or endpoint visualization.



NEB's WarmStart LAMP Kit (DNA & RNA) offers speed and robust sensitivity

A RNA target (HMBS2) was amplified from Jurkat total RNA using the WarmStart LAMP Kit and OptiGene Master Mix (ISO-001). Reactions were performed at 65°C for 74 minutes on a real-time thermocycler (Bio-Rad® CFX96) in triplicate. Time to result is set as the time at which the fluorescence crossed a threshold of 10% of maximal fluorescence. NEB's WarmStart LAMP Kit resulted in faster and more sensitive detection as compared to the OptiGene Master Mix.




optimization tips for LAMP

- Use LAMP primer design software (e.g., Primer Explorer – primerexplorer.jp/e/). Select 2–3 sets for each target and compare performance in a LAMP assay.
- Include loop primers for faster reactions
- Use high magnesium (6–8 mM) and dNTP (1–1.4 mM) concentrations for best results
- Omit betaine, unless it has a demonstrated benefit
- Optimize the reaction temperature (60–65°C for Bst LF and 63–70°C for Bst 2.0/3.0)
- To prevent contamination, use Bst 3.0 or Antarctic Thermo-labile UDG (NEB #M0372), which denatures rapidly

ORDERING INFORMATION:

PRODUCT	NEB #	SIZE
WarmStart LAMP KIT (DNA & RNA)	E1700S/L	1,600/8,000 units
WarmStart Colorimetric LAMP 2X Master Mix (DNA & RNA)	M1800S/L	100/500 rxns
Bst 2.0 WarmStart DNA Polymerase	M0538S/M/L	1,600/8,000 units
WarmStart RTx Reverse Transcriptase	M0380S/L	50/250 rxns



For more detailed information, request the free brochure from your local distributor.


Did you know?


All of the wastewater from our Ipswich, MA campus is treated by a state-of-the-art Solar Aquatics® System. This system utilizes and accelerates the natural purification processes found in streams and wetlands, purifying the water for reuse or for groundwater recharge.


Learn more and take a tour of the facility at www.neb.com/environmentalphilosophy

Tools to support your CRISPR/Cas9

New England Biolabs provides reagents to support a broad variety of CRISPR/Cas9 genome editing approaches. From introduction of Cas9 and single guide RNA (sgRNA) on plasmids, to direct introduction of Cas9 ribonucleoprotein (RNP) and detection of edits using next generation sequencing or enzymatic mutation detection, NEB provides reagents that simplify and shorten genome editing workflows. Generating RNPs for direct introduction requires Cas9 protein and either sgRNA or separate crRNA and tracrRNA. EnGen® Cas9 NLS, *S. pyogenes* is engineered for high genome editing efficiency. The EnGen sgRNA Synthesis Kit combines template assembly and *in vitro* transcription for rapid generation of microgram quantities of custom sgRNA, requiring only a user-supplied single ssDNA oligonucleotide. To determine editing efficiency, the EnGen Mutation Detection Kit provides a full workflow from PCR amplification to T7 Endonuclease I-based mutation detection. Alternatively, NEB supplies Cas9 wild type and restriction enzymes, both of which can be used *in vitro* to determine the extent of editing.


 Schematic representation of Cas9 Nuclease, *S. pyogenes* sequence recognition and DNA cleavage




 Visit www.neb.com/GenomeEditing to view our online tutorial which walks through genome editing with Cas9

EnGen sgRNA Synthesis Kit

The EnGen sgRNA Synthesis Kit simplifies the generation of microgram quantities of custom sgRNAs in an hour or less by combining template synthesis and transcription. The single-tube reaction is easy to set up and requires a single ~55 nt ssDNA target-specific oligonucleotide, which is combined with the Reaction Mix and Enzyme Mix included in the kit. sgRNAs are suitable for use in downstream applications, including CRISPR/Cas9-based genome editing and *in vitro* DNA cleavage. This single-reaction format offers ease-of-use and eliminates separate DNA amplification and template clean up steps.

 EnGen sgRNA Synthesis Kit overview



“This kit is really easy to use and will save us plenty of time in making sgRNAs! Thanks for the streamlined method!”

— Postdoctoral Researcher,
Harvard University


Genome Editing Workflows

EnGen Mutation Detection Kit

A widely used method to identify mutations is the T7 Endonuclease I mutation detection assay. This assay detects heteroduplex DNA that results from the annealing of a DNA strand, including desired mutations, with a wild-type DNA strand. The EnGen Mutation Detection Kit provides optimized reagents for performing robust T7 Endonuclease-based detection of genome editing events.

Workflow for EnGen Mutation Detection Kit

Genomic DNA is amplified with primers bracketing the modified locus. PCR products are then denatured and re-annealed yielding three classes of possible structures. Duplexes containing a mismatch greater than one base are digested by T7 Endonuclease I. The DNA is then electrophoretically separated and fragment analysis is used to estimate targeting efficiency.



online resources

Plasmid Repositories

- addgene.org

CRISPR-gRNA Design Tools

- deskgen.com
- crispr.mit.edu
- zifit.partners.org/ZiFiT
- e-crisp.org
- chopchop.rc.fas.harvard.edu
- benchling.com

Online Forums


- groups.google.com/forum/#!forum/crispr


Organism-specific Resources

- wormcas9hr.weebly.com
- flyrnai.org

Visit www.neb.com/GenomeEditing to find:

- 1 — **Protocols** for applications
- 2 — such as sgRNA synthesis and direct introduction of sgRNA/Cas9 complexes
- 3 —

 **EnGen sgRNA Template Oligo Designer** (accessible through NEBioCalculator® at NEBiocalculator.neb.com)

 Our latest **Genome Editing Brochure** available for download

ORDERING INFORMATION:

PRODUCT	NEB #	SIZE
EnGen sgRNA Synthesis Kit, <i>S. pyogenes</i>	E3322S	20 rxns
EnGen Mutation Detection Kit	E3321S	25 rxns
T7 Endonuclease I	M0302S/L	250/1,250 units
Q5® Hot Start High-Fidelity 2X Master Mix	M0494S/L	100/500 rxns

EnGen Cas9 NLS, <i>S. pyogenes</i>	M0646T/M	400/2,000 pmol
------------------------------------	----------	----------------

***Offer closes June 30th, 2017.**
Please contact your local distributor for details.

SPECIAL PRICE*



For more detailed information, request the free brochure from your local distributor.

LUNA[®]

Universal qPCR & RT-qPCR Reagents



Lighting the way.

Fluorescence-based quantitative real-time qPCR (qPCR) is the "gold standard" for the detection and quantitation of nucleic acids due to its sensitivity and specificity. Luna products from NEB are optimized for qPCR or RT-qPCR, and are available for either intercalating dye or probe-based detection methods. All Luna products provide robust performance on diverse sample sources and target types.

Each Hot Start Taq-based Luna qPCR master mix has been formulated with a unique passive reference dye that is compatible across a wide variety of instrument platforms, including those that require a ROX reference signal. This means that no additional components are required to ensure machine compatibility. The mixes also contain dUTP, enabling carryover prevention when reactions are treated with NEB's Antarctic Thermolabile UDG (NEB #M0372).

A blue visible dye assists in tracking the reagents when pipetting into clear, multi-welled PCR plates.

These features, combined with rapid, sensitive and precise real-time qPCR performance, make Luna the universal choice for all your qPCR experiments.



EXPERIENCE BEST-IN-CLASS PERFORMANCE

- All Luna products have undergone rigorous testing to optimize specificity, sensitivity, accuracy and reproducibility.
- Products perform consistently across a wide variety of sample sources.
- A comprehensive evaluation of commercially-available qPCR and RT-qPCR reagents demonstrates superior performance of Luna products.

OPTIMIZE YOUR RT-qPCR WITH LUNA WARM-START[®] REVERSE TRANSCRIPTASE

- Novel, thermostable reverse transcriptase (RT) improves performance.
- WarmStart RT paired with Hot Start Taq increases reaction specificity and robustness.

MAKE A SIMPLER CHOICE

- One product per application simplifies selection.
- Convenient master mix formats and user-friendly protocols simplify reaction setup.
- Non-interfering, visible tracking dye helps to eliminate pipetting errors.



Download the exclusive article "Development of a High-Throughput Data Analysis Method for Quantitative Real-Time PCR (qPCR)" for additional details and exciting background information on qPCR:

www.neb-online.eu/qPCR-article

Find the right Luna product for your application

		2 Select your detection method	
		Dye-based	Probe-based
1 Select your target	Genomic DNA or cDNA	Luna Universal qPCR Master Mix (NEB #M3003)	Luna Universal Probe qPCR Master Mix (NEB #M3004)
	RNA	Luna Universal One-Step RT-qPCR Kit (NEB #E3005)	Luna Universal Probe One-Step RT-qPCR Kit (NEB #E3006)



ORDERING INFORMATION		
PRODUCT	NEB #	SIZE
Luna [®] Universal qPCR Master Mix	M3003S/L	200/500 rxns
Luna [®] Universal Probe Master Mix	M3004S/L	200/500 rxns
Luna [®] Universal One-Step RT-qPCR Kit	E3005S/L	200/500 rxns
Luna [®] Universal Probe One-Step RT-qPCR Kit	E3006S/L	200/500 rxns

Experience best-in-class performance for your qPCR & RT-qPCR

All NEB products undergo rigorous testing to ensure optimal performance, and Luna is no exception. We took into consideration numerous important traits when evaluating qPCR, including specificity, sensitivity, accuracy and reproducibility, to develop best-in-class qPCR reagents. Furthermore, we did a


comprehensive evaluation of commercially-available qPCR and RT-qPCR reagents, and developed a method of analysis that allows you to quickly compare and evaluate the performance of these products. We wanted to be sure that Luna products will perform to your expectations for all your targets.

Luna products offer exceptional sensitivity, reproducibility and qPCR performance



qPCR targeting human GAPDH was performed using the Luna Universal Probe qPCR Master Mix over a 6-log range of input template concentrations (20 ng – 0.2 pg Jurkat-derived cDNA) with 8 replicates at each concentration. cDNA was generated from Jurkat total RNA using the NEB Proscript® II First Strand cDNA Synthesis Kit (NEB #E6560). NTC = non-template control

Evaluation of commercially-available dye-based qPCR reagents demonstrates the robustness and specificity of Luna



qPCR reagents from NEB and other manufacturers were tested across 16–18 qPCR targets varying in abundance, length and %GC, using either Jurkat genomic DNA or Jurkat-derived cDNA as input (10 genomic DNA targets and 8 cDNA targets on a Bio-Rad real-time instrument, 9 genomic and 7 cDNA targets on an ABI instrument). For each testing condition, data was collected by 2 users and according to manufacturer's specifications. Results were evaluated for efficiency, low input detection and lack of non-template amplification (where ΔC_q = average C_q of lowest input – average C_q of non-template control). In addition, consistency, reproducibility and overall curve quality were assessed (Quality Score). Bar graph indicates % of targets that met acceptable performance criteria (indicated by green box on dot plot and Quality Score > 3). Results for NEB and other major manufacturers are shown: Bio-Rad, SsoAdvanced™ Universal SYBR® Green Supermix; Roche, FastStart™ SYBR Green Master; QIAGEN, QuantiTect® SYBR Green PCR Kit; ABI, PowerUP™ SYBR Green Master Mix; Promega, GoTaq® qPCR Master Mix. NEB's Luna Universal qPCR Master Mix outperformed all other reagents tested.

Visit www.LUNAqPCR.eu for further information, free samples, videos & brochure request.

Request free LUNA® samples on next page →

GET A
**FREE
T-SHIRT***

Request your free **LUNA®**
qPCR or RT-qPCR sample
from your local distributor,
share your feedback and
get a free **LUNA T-Shirt!***



* Please visit www.LUNAqPCR.eu for details. Samples are limited to 2 per work group or as long as stock lasts.

The new Catalog &
Technical Reference 2017/18

Request your free personal copy from
your local distributor and update your
address details!

**COMING
SOON!**

Catalog &
Technical Reference
2017/18

Your local NEB distributor:

BeNeLux:
BIOKÉ
Tel: (+31) 71 720 0220
(BE: 0800 - 71640)
Fax: (+31) 71 891 0019
info@bioke.com
www.bioke.com

BULGARIA
ELTA 90M Ltd.
Tel: (+359) 2 983 9649
Fax: (+359) 2 983 2211
office@elta90.eu
www.elta90.com

CROATIA:
DIAGNOSTICA SKALPELI D.O.O.
Free call: 800-315911
Tel: (+385) 1 377 8484
Fax: (+385) 1 377 8585
skalpeli@skalpeli.hr
www.skalpeli.hr

CZECH REPUBLIC:
BIOTECH A.S.
Tel: (Toll Free) 0800 124683
Fax: 02 72701742
info@ibiotech.cz
www.ibiotech.cz

DENMARK:
BioNordika Denmark A/S
DENMARK
Tel: (39) 56 20 00
Fax: (39) 56 19 42
info@bionordika.dk
www.bionordika.dk

FINLAND:
BioNordika OY
Tel: (+358) 207 410 270
info@bionordika.fi
www.bionordika.fi

GREECE:
BIOLINE SCIENTIFIC
Tel: (0030) 210 522 6547
Fax: (0030) 210 524 4744
info@bioline.gr
www.bioline.gr

HUNGARY:
KVALITEX SCIENTIFIC
Phone: (1) 340 4700
Fax: (1) 339 8274
info@kvalitex.hu
www.kvalitex.hu

ITALY:
EUROCLONE S.P.A.
Free call: 800-315911
Tel: (02) 381951
Fax: (02) 38101465
info@euroclone.it
www.euroclonegroup.it/celbio/

ISRAEL:
ORNAT LTD.
Tel: (+972) 8 947 7077
Fax: (+972) 8 936 3934
ornatbio@ornat.co.il
www.ornat.co.il

NORWAY:
BioNordika AS
Tel +47 23 03 58 00
Fax +47 23 03 58 01
info@bionordika.no
www.bionordika.no

POLAND:
Lab-JOT Ltd. Sp.z o.o. Sp.k.
Stacjonarny: +48 22 335 988 4
Komórkowy: +48 606 338 879
Fax: +48 22 335 981 9
labjot@labjot.com
www.labjot.com

PORTUGAL:
Werfen
Tel: +351 214 247 312
Fax +351 21 425 52 80
dac-pt@werfen.com
http://pt.werfen.com

ROMANIA:
ELTA 90 MEDICAL RESEARCH SRL
Tel: +40 21 2322694
Fax: +40 21 2322696
office@elta90mr.ro
www.elta90mr.ro

RUSSIA:
SkyGen LLC
Tel: +7 495 215 02 22
info@skygen.com
www.skygen.com

SERBIA:
ALFA GENETICS D.O.O.
Tel: 00381 11 2084252
Fax: 00381 11 2084253
office@alfagenetics.rs
www.alfagenetics.rs

SPAIN:
Werfen
Tel: +34 902 20 30 90
Fax: +34 902 22 33 66
customerservice-es@werfen.com
http://es.werfen.com

SLOVENIA:
Mikro+Polo d.o.o.
Tel: +386 2 614 33 00
Fax: +386 2 614 33 20
info@mikro-polo.si
www.mikro-polo.com

SOUTH AFRICA:
INQABA BIOTEC
Tel: +27 123 43 58 29
Fax: +27 86 677 8409
info@inqababiotec.co.za
www.inqababiotec.co.za

SWEDEN, ESTONIA, LATVIA,
LITHUANIA:
BIONORDIKA SWEDEN AB
Tel: +46 8 30 60 10
Fax: +46 8 30 60 15
info@bionordika.se
www.bionordika.se

SWITZERLAND:
BIOCONCEPT
CH 4123 Allschwil 1
Tel: +41 (0)61 486 80 80
Fax: +41 (0)61 486 80 00
info@bioconcept.ch
www.bioconcept.ch

TURKEY:
SACEM HAYAT TEKNOLOJILERI
Gebze Plastikçiler Organize Sanayi
Bölgesi (GEPOSB)
Tel: 0262 751 02 74
Fax: 0262 751 02 75
sacem@sacem.com.tr
www.sacem.com.tr

neb.com

