



BIOKÉ
sharing knowledge

Plesmanlaan 1d
2333 BZ Leiden
The Netherlands
T. +31 (0)71 568 10 00
T. Belgium: 0800 71640
F. +31 (0)71 568 10 10
info@bioke.com
www.bioke.com

PCR Clean-up

User Manual

NucleoSpin[®] 8 Extract II

April 2008 / Rev. 01

MACHERY-NAGEL



Table of contents

1	Components	4
1.1	Kit contents	4
1.2	Reagents and equipment to be supplied by user	4
2	Product description	5
2.1	Basic principle	5
2.2	Kit specifications	5
2.3	General comments	6
2.4	Setting up vacuum	8
2.5	Suitable vacuum manifolds	8
3	Storage conditions and preparation of working solutions	9
4	Safety instructions – risk and safety phrases	10
5	Protocol	11
5.1	General procedure	11
5.2	Standard protocol for DNA clean-up	13
6	Appendix	16
6.1	Troubleshooting	16
6.2	Ordering information	17
6.3	References	17
6.4	Product use restriction / warranty	17

1 Components

1.1 Kit contents

Cat. No.	NucleoSpin® 8 Extract II	
	12 x 8 preps 740668	60 x 8 preps 740668.5
Binding Buffer NT	30 ml	2 x 75 ml
Wash Buffer NT3 (Concentrate) ¹	100 ml	200 ml
Elution Buffer NE ²	25 ml	125 ml
NucleoSpin® Extract Binding Strips (yellow rings)	12	60
MN Wash Plates (including 6 Paper Sheets)	1	5
Tube Strips (with Cap Strips) ³	1	5
Elution Plates, U-bottom (including one Self-adhering PE Foil)	1	5
User Manual	1	1

1.2 Reagents and equipment to be supplied by user

Reagents: 96-100% Ethanol

The **NucleoSpin® 8 Extract II** kit can be used manually with the NucleoVac 96 Vacuum Manifold (Cat. No. 740681) by using the Starter Set A containing Column Holder A and NucleoSpin® Dummy Strips (see ordering information, section 6.1).

¹ For preparation of working solutions and storage conditions see section 3.

² Composition of Elution Buffer NE: 5 mM Tris/HCl, pH 8.5

³ Set of 1 rack with 12 8-well Tube Strips and 12 Cap Strips

2 Product description

2.1 Basic principle

The **NucleoSpin® 8 Extract II** kit is designed for vacuum based clean-up of DNA after amplification in PCR reactions or after enzymatic treatment. Within the **NucleoSpin® 8 Extract II** procedure the addition of chaotropic salt (Buffer NT) leads to a reversible adsorption of the DNA to the silica membrane of the **NucleoSpin® Extract Binding Strips**. Primers, salts, nucleotides, and proteins (polymerases, BSA) are removed in subsequent washing steps using Buffer NT3. Finally, highly pure DNA is eluted in Elution Buffer NE (5 mM Tris/HCl, pH 8.5) or water (pH 8.5), and can be used directly for further applications.

2.2 Kit specifications

- **NucleoSpin® 8 Extract II** is designed for DNA clean-up e.g. purification of PCR products in the convenient 8-well strip format. The 8-well strip format allows for highest flexibility in sample throughput.
- The kit is designed for vacuum use with the NucleoVac 96 Vacuum Manifold and Starter Set A (see ordering information, section 6.2).
- This kit provides reagents and consumables for purification of up to 15 µg highly pure PCR products.
- DNA recovery of 75-90% is obtained for DNA fragments of 64-10,000 bp. Primers, primer-dimers, nucleotides, salts, and polymerase are removed effectively.
- The final concentration of the eluted DNA is 50-150 ng/µl, depending on elution buffer volume used.
- Typically, the A_{260}/A_{280} ratio is >1.8.
- Eluted PCR products are ready-to-use for e.g. automated fluorescent sequencing, cloning or microarray technology.
- Using the **NucleoSpin® 8 Extract II** kit allows simultaneous processing of up to 48 samples typically within 30 minutes.

Kit specification at a glance	
Parameters	NucleoSpin® 8 Extract II
Direct purification of amplified DNA	++
Desalination, removal of enzymes, nucleotides and /or labeling reagents like biotin or radioactive ATP etc.	++
Elution volume	75-175 µl
DNA binding capacity	15 µg
Time/prep	30 min/48 preps

2.3 General comments

- **Membrane capacity:** The silica membrane allows purification of up to 15 µg DNA. The adsorption of DNA to the membrane is pH dependent. High DNA yields are achieved for reaction mixtures with pH <8.0.
- **Elution of purified PCR products:** The efficiency of the DNA elution depends on the pH of the elution buffer. Elution is most effective at pH 8.0-8.5. When using nuclease-free water for elution, the pH value should be checked and, if necessary, adjusted to 8.0-8.5. A lower pH-value of the selected elution buffer may lead to lower recoveries. Yield of larger DNA fragments (>5-10 kbp) can be increased by using prewarmed (70°C) elution buffer. An elution volume of 75-125 µl Buffer NE, as well as a 1-3 min incubation at room temperature of the elution buffer on the silica membrane are recommended.

See table for a correlation between dispensed elution buffer volume and typical recoveries following the standard protocol.

Dispensed elution buffer	75 µl	100 µl	125 µl	150 µl	175 µl
Recovered elution buffer containing PCR-products	30± 5µl	55±5 µl	80±5 µl	105±5 µl	130±5 µl




Fig. 1: Recovery rate and concentration depend on elution volume. Two different PCR products (252 bp, 787 bp) have been purified with NucleoSpin® 8 Extract II kit.

Table 1: Average DNA recovery rate depends on the size of PCR product	
Size of PCR product [bp]	Average DNA recovery rate [%]
64	60 – 80
164	70 – 85
200	70 – 85
490	85 – 95
982	85 – 95
1500	80
2000	75
4000	50 - 60

2.4 Setting up vacuum

Use suitable house vacuum, a membrane pump or a water aspirator. Adjust vacuum as indicated in the individual steps of the protocol. See ordering information, section 6.2, for vacuum regulator.

- Establish a reliable vacuum source for the NucleoVac 96 Vacuum Manifold. The manifold may be used with vacuum pump, house vacuum or water aspirator. We recommend a vacuum of -0.2 to -0.4 bar (pressure difference). Alternatively adjust vacuum that during the purification the sample flows through the column with a rate of 1-2 drops per second.

2.5 Suitable vacuum manifolds

Besides the NucleoVac 96 Vacuum Manifold the NucleoSpin[®] 8 Extract kit can be used with other common vacuum manifolds. For further details see list below:

Vacuum manifold	Compatibility	Additional equipment
Qiagen/ QIAvac 96*	yes	MN Frame (see ordering information)
Promega Vac-Man 96**	yes	-

As for the NucleoVac 96 Vacuum Manifold the Starter Set A is required if other vacuum manifolds are used, too.

*In general the QIAvac 96 is suitable for the use with the Column Holder A and NucleoSpin[®] Extract 8 Binding Strips. Nevertheless, it is recommended to use the MN Frame to adjust the proper height of the MN Wash Plate and Elution Plate (U-bottom), in order to ensure best performance.

**MN Wash Plate can not be used

3 Storage conditions and preparation of working solutions

Attention

Buffer NT contains chaotropic salt! Wear gloves and goggles!

- All components of the **NucleoSpin® 8 Extract II** kit should be stored at room temperature and are stable up to one year.

Before starting any **NucleoSpin® 8 Extract II** protocol prepare the following:

- Add the indicated volume of 96-100 % ethanol to Buffer NT3 Concentrate.

NucleoSpin® 8 Extract II		
Cat. No.	12 x 8 preps	60 x 8 preps
	740 668	740 668.5
Buffer NT3 (Concentrate)	100 ml add 400 ml ethanol	200 ml add 800 ml ethanol

4 Safety instructions – risk and safety phrases

The following components of the **NucleoSpin® 8 Extract II** kits contain hazardous contents.

Wear gloves and goggles and follow the safety instructions given in this section

Contents	Hazard Contents	Hazard Symbol	Risk Phrases	Safety Phrases
NT	Guanidine thiocyanate	 Xn*	Harmful by inhalation, in contact with skin and if swallowed	R 20/21/22 S 13

Risk Phrases

R 20/21/22 Harmful by inhalation, in contact with the skin and if swallowed

Safety Phrases




S 13 Keep away from food, drink and animal feedstuffs

* Label not necessary, if quantity below 125 g or ml (according to 67/548/EEC Art. 25, 1999/45/EC Art. 12 and German GefStoffV § 42 and TRGS 200 7.1)

5 Protocol

5.1 General procedure

For detailed information on each step see section 5.2.

1 Adjust volume of the reaction mix to 100 µl using Tris buffer (pH 7.0-7.5), nuclease-free water (pH 7.0-7.5), or Buffer NE	100 µl reaction volume	
2 Prepare the NucleoVac 96 Vacuum Manifold		
3 Dispense Binding Buffer NT to NucleoSpin® Extract Binding strips	200 µl NT	
4 Transfer PCR samples to NucleoSpin® Extract Binding Strips and mix	100 µl diluted sample, mix	
5 Bind DNA to silica membrane of the NucleoSpin® Extract Binding Strips by applying vacuum	-0.2 to -0.4 bar* 1 min	
6 Wash silica membrane	2 x 900 µl NT3	Column Holder A with NucleoSpin® Extract Binding Strips and MN Wash Plate
	-0.2 to -0.4 bar 1 min	
7 Remove MN Wash Plate		
8 Dry silica membrane of the NucleoSpin® Extract Binding Strips	10-15 min, max. vacuum	
<i>Optional: Before applying vacuum dry the outlets of the NucleoSpin® Extract Binding Strips by placing them on a Paper Sheet.</i>		Column Holder A with NucleoSpin® Extract Binding Strips


* Reduction of atmospheric pressure

9 Insert rack with Tube Strips or Elution Plate, U-bottom

10 Elute DNA

75-150 µl NE

-0.2 to -0.4 bar*
1 - 2 min



Column Holder A with
NucleoSpin® Extract
Binding Strips and
Tube Strips or Elution
Plate, U-bottom

* Reduction of atmospheric pressure

5.2 Standard protocol for DNA clean-up

1 Adjust volume of sample reaction mix to 100 µl

For sample volumes of less than 100 µl: Before starting the purification procedure, add Tris buffer (10 mM, pH 7.0), nuclease-free water (pH 7.0-7.5), or Buffer NE to **adjust the reaction mixture to a final volume of 100 µl**.

Note: Adjusting the reaction volume to 100 µl allows easier and more convenient transfer of the samples with a multichannel pipette.

2 Prepare the NucleoVac 96 Vacuum Manifold

Insert an appropriate number of NucleoSpin® Extract Binding Strips (yellow rings) into a Column Holder A. Close any unused openings of the Column Holder A with NucleoSpin® Dummy Strips.

Note: Make sure that the NucleoSpin® Extract Binding Strips are inserted tightly into the column holder. Uneven or not properly inserted strips may prevent sealing when vacuum is applied to the manifold.

Insert spacers (*MTP/Multi-96 Plate*), notched side up, into the grooves located on the short sides of the manifold. Insert deep waste reservoir into the center of the manifold. Place the MN Wash Plate on top on the spacers. Insert Column Holder A with inserted NucleoSpin® Extract Binding Strips into the manifold lid and place lid on the manifold base. Check and adjust the vacuum (pressure difference -0.2 bar). Release the vacuum.

3 Dispense Binding Buffer NT to the NucleoSpin® Extract Binding Strips

Add 2 volumes of Buffer NT for one volume of sample (e.g. 200 µl Buffer NT for 100 µl reaction mix).

4 Transfer samples to the NucleoSpin® Extract Binding Strips and mix

Mix sample and Buffer NT by pipetting up and down 5 times.

Optional: Premix Buffer NT and sample in a suitable 96-well plate (e.g. Round-well Block). Transfer the mixture to the NucleoSpin® Extract Binding Strips.

5 Bind DNA to silica membrane of the NucleoSpin® Extract Binding Strips

Apply vacuum by opening the valve and press down the Column Holder A slightly until flow-through starts. Allow the samples to pass the wells and release the vacuum.

6 Wash silica membrane of the NucleoSpin® Extract Binding Strips

Add **900 µl** of **Buffer NT3** (with ethanol added) to each well of the NucleoSpin® Extract Binding Strips. Apply vacuum and allow the buffer to pass the wells. Release the vacuum.

Repeat this washing step once.

7 Remove MN Wash Plate

After the final washing step close the valve, release the vacuum and remove the Column Holder A with inserted NucleoSpin® Extract Binding Strips. Remove manifold lid, MN Wash Plate, and waste container from the vacuum manifold.

8 Dry silica membrane of the NucleoSpin® Extract Binding Strips

Remove any residual washing buffer from the NucleoSpin® Extract Binding Strips. If necessary, dry the nozzles of the strips on a clean paper sheet (supplied with the kit) or soft tissue until no further drops come out. Insert the Column Holder A holding the NucleoSpin® Extract Binding Strips into the lid and close the manifold. Apply maximum vacuum for 10-15 min to dry the membrane completely. This step is necessary to eliminate traces of ethanol

Note: The ethanol in Buffer NT3 inhibits enzymatic reactions and has to be removed completely before eluting DNA.

Finally, release the vacuum.

9 Insert rack with Tube Strips for elution

Insert spacers (*Microtube Rack*), notched side up, into the grooves located at the short sides of the vacuum manifold. Place the Tube Strip rack on the spacers inside the manifold base and insert the appropriate number of Tube Strips. Insert the column holder with the NucleoSpin® Extract Binding Strips in the manifold lid and place it on the manifold base.

Alternatively, elution can be performed into a 96-well microtiter plate (supplied with the kit). Insert spacers (*MTP/Multi-96 Plate*), notched side up, into the grooves located at the short sides of the vacuum manifold. Place the Elution Plate on the spacers inside the manifold base. Insert the column holder with the NucleoSpin® Extract Binding Strips in the manifold lid and place it on the manifold base.

10 Elute DNA

Elute the DNA by adding 75-150 µl of Buffer NE (5 mM Tris-HCl, pH 8.5) or nuclease free water (pH 8.5) to each well of the NucleoSpin® Extract Binding Strips .

The elution buffer should be dispensed carefully onto the center of the silica membrane. Incubate the buffer on the membrane for 1-3 minutes. Apply vacuum (-0.2 to -0.4 bar). Press down the column holder slightly and collect the eluted DNA. After the elution buffer has passed the wells release the vacuum.

Remove the Rack with Tube Strips or Elution Plate containing the eluted DNA.

6 Appendix

6.1 Troubleshooting

Problem	Possible cause and suggestions
Poor DNA yield	<i>Buffer NT3 did not contain ethanol</i>
	<ul style="list-style-type: none"> • Addition of the indicated volume of 96-100 % ethanol to the Buffer NT3 Concentrate is required before use.
	<i>Inappropriate elution buffer</i>
<ul style="list-style-type: none"> • Elution is most effective at pH 8.5. If using water, check the pH value before use. Always use low salt Buffer NE or water for elution. 	
<i>Elution buffer volume insufficient</i>	<ul style="list-style-type: none"> • Optimal elution is achieved for an elution buffer volume of 100-150 µl. Do not use less than 75 µl of elution buffer.
Suboptimal performance of DNA in sequencing reactions	<i>Carryover of ethanol</i>
	<ul style="list-style-type: none"> • Be sure to remove all of ethanolic Buffer NT3 after the final wash step. Dry NucleoSpin[®] Extract Binding Strips for at least 10 min with maximum airflow.
	<i>Elution of DNA with TE buffer</i>
<ul style="list-style-type: none"> • EDTA may inhibit sequencing reactions. Repurify or precipitate DNA by ethanol and elute/redissolve in Buffer NE or in water. 	
<i>Eluted DNA contains residual primers/primer dimers</i>	<ul style="list-style-type: none"> • Minimize amount of primers in PCR reaction mixture. Make sure that the ratio of binding buffer NT:PCR reaction mixture is 2:1.

6.2 Ordering information

Product	Cat. No.	Pack of
NucleoSpin® 8 Extract II	740668	12 x 8 preps
NucleoSpin® 8 Extract II	740668.5	60 x 8 preps
Round-well Block	740671	20
Elution Plate, U-bottom	740672	20
Self-adhering PE Foil	740676	50 sheets
Tube Strips	740637	5 racks
Cap Strips	740638	30
NucleoVac 96 Vacuum Manifold	740681	1
NucleoVac Vacuum Regulator	740641	1
Starter Set A (for use of the NucleoSpin® 8-well strips on the NucleoVac 96 Vacuum Manifold)	740684	1 set

6.3 References

Vogelstein, B. & Gillespie, D. (1979) Proc. Natl. Acad. Sci. USA 76, 615-619.

6.4 Product use restriction / warranty

NucleoSpin® 8 Extract II kit components were developed, designed, distributed, and sold **for RESEARCH PURPOSES ONLY**. They are suitable **for IN - VITRO USES only**. Furthermore is no claim or representation intended for its use to identify any specific organism or for clinical use (diagnostic, prognostic, therapeutic, or blood banking).

It is rather in the responsibility of the user to verify the use of the **NucleoSpin® 8 Extract II** kits for a specific application range as the performance characteristic of this kit has not been verified to a specific organism.

This MACHEREY-NAGEL product is shipped with documentation stating specifications and other technical information. MACHEREY-NAGEL guarantees to meet the stated specifications. MACHEREY-NAGEL's sole obligation and the customer's sole remedy is limited to replacement of products free of charge in the event products fail to perform as warranted. Supplementary reference is made to the

general business terms and conditions of MACHEREY-NAGEL, which are printed on the price list. Please contact us if you wish an extra copy.

MACHEREY-NAGEL does not warrant against damages or defects arising in shipping and handling (transport insurance for customers excluded), or out of accident or improper or abnormal use of this product; against defects in products or components not manufactured by MACHEREY-NAGEL, or against damages resulting from such non-MACHEREY-NAGEL components or products.

MACHEREY-NAGEL makes no other warranty of any kind whatsoever, and SPECIFICALLY DISCLAIMS AND EXCLUDES ALL OTHER WARRANTIES OF ANY KIND OR NATURE WHATSOEVER, DIRECTLY OR INDIRECTLY, EXPRESS OR IMPLIED, INCLUDING, WITHOUT LIMITATION, AS TO THE SUITABILITY, PRODUCTIVITY, DURABILITY, FITNESS FOR A PARTICULAR PURPOSE OR USE, MERCHANTABILITY, CONDITION, OR ANY OTHER MATTER WITH RESPECT TO MACHEREY-NAGEL PRODUCTS.

In no event shall MACHEREY-NAGEL be liable for claims for any other damages, whether direct, incidental, foreseeable, consequential, or special (including but not limited to loss of use, revenue or profit), whether based upon warranty, contract, tort (including negligence) or strict liability arising in connection with the sale or the failure of MACHEREY-NAGEL products to perform in accordance with the stated specifications.

The warranty provided herein and the data, specifications and descriptions of this MACHEREY-NAGEL product appearing in MACHEREY-NAGEL published catalogues and product literature are MACHEREY-NAGEL's sole representations concerning the product and warranty. No other statements or representations, written or oral, by MACHEREY-NAGEL's employees, agent or representatives, except written statements signed by a duly authorized officer of MACHEREY-NAGEL are authorized; they should not be relied upon by the customer and are not a part of the contract of sale or of this warranty.

Product claims are subject to change. Therefore please contact our Technical Service Team for the most up-to-date information on MACHEREY-NAGEL products. You may also contact your local distributor for general scientific information. Applications mentioned in MACHEREY-NAGEL literature are provided for informational purposes only. MACHEREY-NAGEL does not warrant that all applications have been tested in MACHEREY-NAGEL laboratories using MACHEREY-NAGEL products. MACHEREY-NAGEL does not warrant the correctness of any of those applications.

Please contact:

MACHEREY-NAGEL Germany
Tel.: +49 (0) 2421 969 270
e-mail: TECH-BIO@mn-net.com

Last updated 12/2006, Rev. 02