#14166 Store at RT

Hematoxylin

www.cellsignal.com

Support: 877-678-TECH (8324) www.cellsignal.com/support


Orders: 877-616-CELL (2355) orders@cellsignal.com

500 ml

New 11/14

For Research Use Only. Not For Use In Diagnostic Procedures.

Description: Hematoxylin is a blue nuclear counterstain for use in immunohistochemical assays. It yields crisp staining detail with superior contrast when used in conjunction with SignalStain® DAB Substrate Kit #8059. It is also compatible with SignalStain® Mounting Medium #14177.



A) Solutions and Reagents

Supplied Reagents


1. Hematoxylin #14166

Additional Reagents (Not Supplied)

- 1. Acid Rinse Solution: To prepare 100 ml, add 2 ml of glacial acetic acid to 98 ml of dH₂0.
- 2. Bluing Solution: To prepare 100 ml, add 1.5 ml of 30% ammonium hydroxide (NH $_{\rm A}$ OH) to 98.5ml of 70% ethanol

B) Nuclear Counterstaining

- 1. Rinse slides in tap water.
- 2. Completely submerse slides in Hematoxylin #14166 and incubate for 1-5 minutes.
- 3. Wash slides with tap water until rinse water is clear.
- 4. Immerse slides 10 times in acid rinse solution.
- 5. Rinse slides 10 times in tap water.
- 6. Immerse slides for 1 minute in bluing solution
- 7. Rinse slides 10 times in tap water.
- 8. Mount sections with SignalStain® Mounting Medium #14177.



Immunohistochemical analysis of paraffin-embedded human lung carcinoma using Claudin-1 (D5H1D) XP® Rabbit mAb #13255 (brown) followed by a counterstain with Hematoxylin (blue).

Storage: Supplied as a working solution. Store at room temperature. Product is stable for 12 months upon receipt when stored properly.

Thank you for your recent purchase. If you would like to provide a review visit www.cellsignal.com/comments.

