




NextGENe[®]

Viewer Setup Guide

Local licensed Viewer

1. Log on to the computer as a Local Administrator. In order to meet User Account Control settings, you may also need to right-click on the NextGENe Viewer setup and choose Run As Administrator.
2. Choose to install NextGENe Viewer.



1. Open NextGENe Viewer (as Administrator). A Warning dialog indicates that a valid license cannot be found.
2. For Local Installation, select "Register".
3. Input your Account, Password and Email and click Register.
4. Click OK when prompted with successful registration.

The screenshot shows the NextGENe Viewer application window with a menu bar (File, Process, Paired View, Reports, Search, Tools, Help) and a toolbar. A "Warning" dialog box is open, displaying a yellow warning icon and the following text: "This application has not been registered or the license has expired. Click 'Register' to complete this process. If you have network license and installed License Manager S please click 'Network Configure'". Below the text are three buttons: "Register", "Network Configure", and "Cancel".

Overlaid on the main window is a registration dialog box. It contains the text: "Register NextGENe Viewer using the account and password from email or as printed on the installation disc. For support please contact tech_support@softgenetics.com". Below this text is a link "Offline Registration >>". There are four input fields: "Request Code" (containing "MUVUTXhFRFF3RWRReUWVTzVRVWJDUKRS"), "Account", "Password", and "Email". At the bottom of the dialog are "Register" and "Cancel" buttons.

At the bottom of the screen, an "Information" dialog box is open, displaying an information icon and the text: "You have register your application successfully." Below the text is an "OK" button.