e at 4°C	PathScan <sup>®</sup> Total ELISA Kit						
Store						Orders:	877-616-CELL (2355) orders@cellsignal.com
82						Support:	877-678-TECH (8324)
#14784	1 Kit (96 assays) Species Cross Reactivity	UniProt ID:	Entrez-Gene Id:			Web:	info@cellsignal.com cellsignal.com
#	H	#Q9NZQ7	#29126		3 Trask	Lane   Danvers   M	lassachusetts   01923   USA
For Re	esearch Use Only. Not for U	se in Diagnosti	ic Procedures.				1 I
Produ	ict Includes			Product #	Quantity	Color	Storage Temp
	Rabbit mAb Coated Microwells			2/1721	96 tests		٨°C

PD-L1 Rabbit mAb Coated Microwells	24721	24721 96 tests		4°C	
PD-L1 Mouse Detection mAb	14785	1 ea	Green (Lyophilized)	4°C	
Anti-mouse IgG, HRP-linked Antibody (ELISA Formulated)	13304	1 ea	Red (Lyophilized)	4°C	
Detection Antibody Diluent	13339	11 ml	Green	4°C	
HRP Diluent	13515	11 ml	Red	4°C	
TMB Substrate	7004	11 ml		4°C	
STOP Solution	7002	11 ml		4°C	
Sealing Tape	54503	2 ea		4°C	
ELISA Wash Buffer (20X)	9801	25 ml		4°C	
ELISA Sample Diluent	11083	25 ml	Blue	4°C	
Cell Lysis Buffer (10X)	9803	15 ml		-20°C	

\*The microwell plate is supplied as 12 8-well modules - Each module is designed to break apart for 8 tests.

Description	The PathScan <sup>®</sup> Total PD-L1 Sandwich ELISA Kit is a solid phase sandwich enzyme-linked immunosorbent assay (ELISA) that detects endogenous levels of PD-L1 protein. A rabbit antibody toward intracellular PD-L1 has been coated onto the microwells. After incubation with cell lysates, PD-L1 protein is captured by the coated antibody. Following extensive washing, a mouse detection antibody toward extracellular PD-L1 is added to detect the captured PD-L1 protein. Anti-mouse IgG, HRP-linked antibody is then used to recognize the bound detection antibody. HRP substrate, TMB, is added to develop color. The magnitude of the absorbance for the developed color is proportional to the quantity of PD-L1 protein.
	*Antibodies in this kit are custom formulations specific to kit.
Specificity/Sensitivity	PathScan <sup>®</sup> Total PD-L1 Sandwich ELISA Kit detects endogenous levels of PD-L1 protein as shown in Figure 1. The kit sensitivity is shown in Figure 2. This kit detects proteins from the indicated species, as determined through in-house testing, but may also detect homologous proteins from other species.
Background	Programmed cell death 1 ligand 1 (PD-L1, B7-H1, CD274) is a member of the B7 family of cell surface ligands that regulate T cell activation and immune responses. The PD-L1 ligand binds the PD-1 transmembrane receptor and inhibits T cell activation. PD-L1 was discovered following a search for novel B7 protein homologs and was later shown to be expressed by antigen presenting cells, activated T cells, and tissues including placenta, heart, and lung (1-3). Similar in structure to related B7 family members, PD-L1 protein contains extracellular IgV and IgC domains and a short, cytoplasmic region. Research studies demonstrate that PD-L1 is expressed in several tumor types, including melanoma, ovary, colon, lung, breast, and renal cell carcinomas (4-6). Expression of PD-L1 in cancer is associated with tumor-infiltrating lymphocytes, which mediate PD-L1 expression through the release of interferon gamma (7). Additional research links PD-L1 expression to cancers associated with viral infections (8,9).

1/1/24, 7:00 AM	PathScan® Total PD-L1 Sandwich ELISA Kit (#14784) Datasheet Without Images Cell Signaling Technology
Background References	<ol> <li>Dong, H. et al. (1999) Nat Med 5, 1365-9.</li> <li>Freeman, G.J. et al. (2000) J Exp Med 192, 1027-34.</li> <li>Liang, S.C. et al. (2003) Eur J Immunol 33, 2706-16.</li> <li>Dong, H. et al. (2002) Nat Med 8, 793-800.</li> <li>Thompson, R.H. et al. (2006) Cancer Res 66, 3381-5.</li> <li>Pardoll, D.M. (2012) Nat Rev Cancer 12, 252-64.</li> <li>Taube, J.M. et al. (2012) Sci Transl Med 4, 127ra37.</li> <li>Lyford-Pike, S. et al. (2013) Cancer Res 73, 1733-41.</li> <li>Chen, B.J. et al. (2013) Clin Cancer Res 19, 3462-73.</li> </ol>
Cross-Reactivity	<ul> <li>Key</li> <li>H: human M: mouse R: rat Hm: hamster Mk: monkey Vir: virus Mi: mink C: chicken Dm: D. melanogaster</li> <li>X: Xenopus Z: zebrafish B: bovine Dg: dog Pg: pig Sc: S. cerevisiae Ce: C. elegans Hr: horse</li> <li>GP: Guinea Pig Rab: rabbit All: all species expected</li> </ul>
Trademarks and Patents	Cell Signaling Technology is a trademark of Cell Signaling Technology, Inc. E1L3N is a registered trademark of Cell Signaling Technology, Inc. PathScan is a registered trademark of Cell Signaling Technology, Inc. XP is a registered trademark of Cell Signaling Technology, Inc. All other trademarks are the property of their respective owners. Visit cellsignal.com/trademarks for more information.
Limited Uses	Except as otherwise expressly agreed in a writing signed by a legally authorized representative of CST, the following terms apply to Products provided by CST, its affiliates or its distributors. Any Customer's terms and conditions that are in addition to, or different from, those contained herein, unless separately accepted in writing by a legally authorized representative of CST, are rejected and are of no force or effect.
	Products are labeled with For Research Use Only or a similar labeling statement and have not been approved, cleared, or licensed by the FDA or other regulatory foreign or domestic entity, for any purpose. Customer shall not use any Product for any diagnostic or therapeutic purpose, or otherwise in any manner that conflicts with its labeling statement. Products sold or licensed by CST are provided for Customer as the end-user and solely for research and development uses. Any use of Product for diagnostic, prophylactic or therapeutic purposes, or any purchase of Product for reseale (alone or as a component) or other commercial purpose, requires a separate license from CST. Customer shall (a) not sell, license, loan, donate or otherwise transfer or make available any Product to any third party, whether alone or in combination with other materials, or use the Products to manufacture any commercial products, (b) not copy, modify, reverse engineer, decompile, disassemble or otherwise attempt to discover the underlying structure or technology of the Products, or use the Products for the purpose of developing any products any trademarks, trade names, logos, patent or copyright notices or markings, (d) use the Products solely in accordance with CST Product Terms of Sale and any applicable documentation, and (e) comply with any license, terms of service or similar agreement with respect to any third party products or services used by Customer in connection with the Products.

# #**14784** PathScan<sup>®</sup> Total PD-L1 Sandwich ELISA Kit



# **ELISA Colorimetric (Lyophilized)**

#### A. Solutions and Reagents

NOTE: Prepare solutions with purified water.

- 1. **Microwell strips**: Bring all to room temperature before use.
- 2. Detection Antibody: Supplied lyophilized as a green colored cake or powder. Add 1.0 ml of Detection Antibody Diluent (green solution) to yield a concentrated stock solution. Incubate at room temperature for 5 min with occasional gentle mixing to fully reconstitute. To make the final working solution, add the full 1.0 ml volume of reconstituted Detection Antibody to 10.0 ml of Detection Antibody Diluent in a clean tube and gently mix. Unused working solution may be stored for 4 weeks at 4°C.
- 3. HRP-Linked Antibody\*: Supplied lyophilized as a red colored cake or powder. Add 1.0 ml of HRP Diluent (red solution) to yield a concentrated stock solution. Incubate at room temperature for 5 min with occasional gentle mixing to fully reconstitute. To make the final working solution, add the full 1.0 ml volume of reconstituted HRP-Linked Antibody to 10.0 ml of HRP Diluent in a clean tube and gently mix. Unused working solution may be stored for 4 weeks at 4°C.
- 4. Detection Antibody Diluent: Green colored diluent for reconstitution and dilution of the detection antibody (11 ml provided).
- 5. HRP Diluent: Red colored diluent for reconstitution and dilution of the HRP-Linked Antibody (11 ml provided).
- 6. Sample Diluent: Blue colored diluent provided for dilution of cell lysates.
- 7. 1X Wash Buffer: Prepare by diluting 20X Wash Buffer (included in each PathScan<sup>®</sup> Sandwich ELISA Kit) in purified water.
- 8. **Cell Lysis Buffer**: 10X Cell Lysis Buffer #9803: This buffer can be stored at 4°C for short-term use (1–2 weeks). Recommended: Add 1 mM phenylmethylsulfonyl fluoride (PMSF) immediately before use.
- 9. TMB Substrate (#7004).
- 10. STOP Solution (#7002).

\*NOTE: Some PathScan<sup>®</sup> ELISA Kits may include HRP-Linked Streptavidin in place of HRP-Linked Antibody.

### **B. Preparing Cell Lysates**

#### For adherent cells.

- 1. Aspirate media when the culture reaches 80-90% confluence. Treat cells by adding fresh media containing regulator for desired time.
- 2. Remove media and rinse cells once with ice-cold 1X PBS.
- 3. Remove PBS and add 0.5 ml ice-cold 1X Cell Lysis Buffer plus 1 mM PMSF to each plate (10 cm diameter) and incubate the plate on ice for 5 min.
- 4. Scrape cells off the plate and transfer to an appropriate tube. Keep on ice.
- 5. Sonicate lysates on ice.
- 6. Microcentrifuge for 10 min (x14,000 rpm) at 4°C and transfer the supernatant to a new tube. The supernatant is the cell lysate. Store at -80°C in single-use aliquots.

#### For suspension cells

- 1. Remove media by low speed centrifugation (~1200 rpm) when the culture reaches 0.5–1.0 x 10<sup>6</sup> viable cells/ml. Treat cells by adding fresh media containing regulator for desired time.
- 2. Collect cells by low speed centrifugation (~1200 rpm) and wash once with 5–10 ml ice-cold 1X PBS.
- 3. Cells harvested from 50 ml of growth media can be lysed in 2.0 ml of 1X Cell Lysis Buffer plus 1 mM PMSF.
- 4. Sonicate lysates on ice.
- 5. Microcentrifuge for 10 min (x14,000 rpm) at 4°C and transfer the supernatant to a new tube. The supernatant is the cell lysate. Store at -80°C in single-use aliquots.

## C. Test Procedure

- 1. After the microwell strips have reached room temperature, break off the required number of microwells. Place the microwells in the strip holder. Unused microwells must be resealed and stored at 4°C immediately.
- 2. Cell lysates can be undiluted or diluted with Sample Diluent (supplied in each PathScan<sup>®</sup> Sandwich ELISA Kit, blue color). Individual datasheets for each kit provide a sensitivity curve that serves as a reference for selection of an appropriate starting lysate
- concentration. The sensitivity curve shows typical kit assay results across a range of lysate concentration points. 3. Add 100  $\mu$ l of each undiluted or diluted cell lysate to the appropriate well. Seal with tape and press firmly onto top of microwells.
- Incubate the plate for 2 hr at 37°C. Alternatively, the plate can be incubated overnight at 4°C.
- 4. Gently remove the tape and wash wells:
  - 1. Discard plate contents into a receptacle.
  - 2. Wash 4 times with 1X Wash Buffer, 200 µl each time for each well.
  - 3. For each wash, strike plates on fresh towels hard enough to remove the residual solution in each well, but do not allow wells to completely dry at any time.
  - 4. Clean the underside of all wells with a lint-free tissue.
- 5. Add 100 μl of reconstituted Detection Antibody (green color) to each well (refer to Section A, Step 2). Seal with tape and incubate the plate at 37°C for 1 hr.

#### 1/1/24, 7:00 AM PathScan® Total PD-L1 Sandwich ELISA Kit (#14784) Datasheet Without Images Cell Signaling Technology

- 6. Repeat wash procedure (Section C, Step 4).
- 7. Add 100 µl of reconstituted HRP-Linked secondary antibody (red color) to each well (refer to Section A, Step 3). Seal with tape and incubate the plate for 30 min at 37°C.
- 8. Repeat wash procedure (Section C, Step 4).
- 9. Add 100 µl of TMB Substrate to each well. Seal with tape and incubate the plate for 10 min at 37°C or 30 min at 25°C.
- 10. Add 100 µl of STOP Solution to each well. Shake gently for a few seconds.

NOTE: Initial color of positive reaction is blue, which changes to yellow upon addition of STOP Solution.

- 11. Read results.

  - Visual Determination: Read within 30 min after adding STOP Solution.
     Spectrophotometric Determination: Wipe underside of wells with a lint-free tissue. Read absorbance at 450 nm within 30 min after adding STOP Solution.

posted November 2013

Orders: 877-616-CELL (2355) • orders@cellsignal.com • Support: 877-678-TECH (8324) • info@cellsignal.com • Web: cellsignal.com For Research Use Only. Not for Use in Diagnostic Procedures.