

#89947 Store at -20°C

ER Stress-induced Autophagy Antibody Sampler Kit

1 Kit (7 x 20 microliters)



Cell Signaling
TECHNOLOGY®

Orders: 877-616-CELL (2355)
orders@cellsignal.com

Support: 877-678-TECH (8324)

Web: info@cellsignal.com
cellsignal.com

3 Trask Lane | Danvers | Massachusetts | 01923 | USA

For Research Use Only. Not for Use in Diagnostic Procedures.

Product Includes	Product #	Quantity	Mol. Wt	Isotype/Source
BiP (C50B12) Rabbit mAb	3177	20 µl	78 kDa	Rabbit IgG
eIF2α (D7D3) XP® Rabbit mAb	5324	20 µl	38 kDa	Rabbit IgG
Phospho-eIF2α (Ser51) (D9G8) XP® Rabbit mAb	3398	20 µl	38 kDa	Rabbit IgG
Atg12 (D88H11) Rabbit mAb	4180	20 µl	16, 55 kDa	Rabbit IgG
Beclin-1 (D40C5) Rabbit mAb	3495	20 µl	60 kDa	Rabbit IgG
JNK1 (2C6) Mouse mAb	3708	20 µl	46, 54 kDa	Mouse IgG1
Phospho-SAPK/JNK (Thr183/Tyr185) (81E11) Rabbit mAb	4668	20 µl	46, 54 kDa	Rabbit IgG
Anti-rabbit IgG, HRP-linked Antibody	7074	100 µl		Goat
Anti-mouse IgG, HRP-linked Antibody	7076	100 µl		Horse

Please visit [cellsignal.com](https://www.cellsignal.com) for individual component applications, species cross-reactivity, dilutions, protocols, and additional product information.

Description

The ER Stress-induced Antibody Sampler Kit contains reagents to investigate ER stress-induced signaling within the cell. The kit contains enough primary antibodies to perform four western blot experiments per primary antibody.

Storage

Supplied in 10 mM sodium HEPES (pH 7.5), 150 mM NaCl, 100 µg/ml BSA, 50% glycerol and less than 0.02% sodium azide. Store at -20°C. Do not aliquot the antibodies.

Background

The endoplasmic reticulum (ER) is an organelle with essential biosynthetic and signaling functions in eukaryotic cells (1). Post synthesis of secretory and transmembrane proteins on polysomes, proteins are translocated into the ER where they are often modified by disulfide bond formation, amino-linked glycosylation, and folding. Different physiological and pathological conditions can disturb proper protein folding in the ER causing ER stress (1). ER stress activates an intracellular signaling transduction pathway called unfolded protein response (UPR) and autophagy to avoid cell death (2). The main role of UPR is to improve the protein load on the ER by shutting down protein translation and gene transcription to enhance ER's folding capacity (2). On the other hand, autophagy is a catabolic process for the autophagosomal-lysosomal degradation of bulk cytoplasmic contents (3,4). One of the chaperones aiding in proper protein folding is Binding immunoglobulin Protein (BiP) (5,6). BiP works by binding to misfolded proteins to prevent them from forming aggregates and assists in proper refolding (7). The molecular machinery of autophagy was largely discovered in yeast and referred to as autophagy-related (Atg) genes. Formation of the autophagosome involves a ubiquitin-like conjugation system in which Atg12 is covalently bound to Atg5 and targeted to autophagosome vesicles (8-10). One of the proteins critical to autophagy process is Beclin-1, the mammalian orthologue of the yeast autophagy protein Apg6/Vps30 (11). Beclin-1 can complement defects in yeast autophagy caused by loss of Apg6 and can also stimulate autophagy when overexpressed in mammalian cells (12). Mammalian Beclin-1 was originally isolated in a yeast two-hybrid screen for Bcl-2 interacting proteins and has been shown to interact with Bcl-2 and Bcl-xL, but not with Bax or Bak (13). Phosphorylation of the eukaryotic initiation factor 2 (eIF2) α subunit is a well-documented mechanism to downregulate protein synthesis under a variety of stress conditions. eIF2 binds GTP and Met-tRNAi and transfers Met-tRNA to the 40S subunit to form the 43S preinitiation complex (14,15). Kinases that are activated by viral infection (PKR) can phosphorylate the α subunit of eIF2 (16,17). Induction of PKR by IFN-γ and TNF-α induces potent phosphorylation of eIF2α at Ser51 (18,19). There are three SAPK/JNK genes each of which undergoes alternative splicing, resulting in numerous isoforms (20). The IRE1, a transmembrane serine/threonine kinase (21,22), through its kinase activity activates SAPK/JNK in the early stage of ER stress in order to induce autophagosome formation (23).

Background References

1. Verfaillie, T. et al. (2010) *Int J Cell Biol* 2010, 930509.
2. Ogata, M. et al. (2006) *Mol Cell Biol* 26, 9220-31.
3. Reggiori, F. and Klionsky, D.J. (2002) *Eukaryot Cell* 1, 11-21.
4. Codogno, P. and Meijer, A.J. (2005) *Cell Death Differ* 12 Suppl 2, 1509-18.
5. Wabl, M. and Steinberg, C. (1982) *Proc Natl Acad Sci U S A* 79, 6976-8.

6. Haas, I.G. and Wabl, M. (2002) *Nature* 306, 387-9.
7. Kohno, K. et al. (1993) *Mol Cell Biol* 13, 877-90.
8. Mizushima, N. et al. (1998) *J Biol Chem* 273, 33889-92.
9. Mizushima, N. et al. (1998) *Nature* 395, 395-8.
10. Suzuki, K. et al. (2001) *EMBO J* 20, 5971-81.
11. Kametaka, S. et al. (1998) *J Biol Chem* 273, 22284-91.
12. Liang, X.H. et al. (1999) *Nature* 402, 672-6.
13. Liang, X.H. et al. (1998) *J Virol* 72, 8586-96.
14. Kimball, S.R. (1999) *Int J Biochem Cell Biol* 31, 25-9.
15. de Haro, C. et al. (1996) *FASEB J* 10, 1378-87.
16. Kaufman, R.J. (1999) *Genes Dev* 13, 1211-33.
17. Sheikh, M.S. and Fornace, A.J. (1999) *Oncogene* 18, 6121-8.
18. Cheshire, J.L. et al. (1999) *J Biol Chem* 274, 4801-6.
19. Zamanian-Daryoush, M. et al. (2000) *Mol Cell Biol* 20, 1278-90.
20. Kyriakis, J.M. and Avruch, J. (2001) *Physiol Rev* 81, 807-69.
21. Nikawa, J. and Yamashita, S. (1992) *Mol Microbiol* 6, 1441-6.
22. Cox, J.S. et al. (1993) *Cell* 73, 1197-206.
23. Urano, F. et al. (2000) *Science* 287, 664-6.

Trademarks and Patents

Cell Signaling Technology is a trademark of Cell Signaling Technology, Inc.

XP is a registered trademark of Cell Signaling Technology, Inc.

U.S. Patent No. 7,429,487, foreign equivalents, and child patents deriving therefrom.

All other trademarks are the property of their respective owners. Visit cellsignal.com/trademarks for more information.

Limited Uses

Except as otherwise expressly agreed in a writing signed by a legally authorized representative of CST, the following terms apply to Products provided by CST, its affiliates or its distributors. Any Customer's terms and conditions that are in addition to, or different from, those contained herein, unless separately accepted in writing by a legally authorized representative of CST, are rejected and are of no force or effect.

Products are labeled with For Research Use Only or a similar labeling statement and have not been approved, cleared, or licensed by the FDA or other regulatory foreign or domestic entity, for any purpose. Customer shall not use any Product for any diagnostic or therapeutic purpose, or otherwise in any manner that conflicts with its labeling statement. Products sold or licensed by CST are provided for Customer as the end-user and solely for research and development uses. Any use of Product for diagnostic, prophylactic or therapeutic purposes, or any purchase of Product for resale (alone or as a component) or other commercial purpose, requires a separate license from CST. Customer shall (a) not sell, license, loan, donate or otherwise transfer or make available any Product to any third party, whether alone or in combination with other materials, or use the Products to manufacture any commercial products, (b) not copy, modify, reverse engineer, decompile, disassemble or otherwise attempt to discover the underlying structure or technology of the Products, or use the Products for the purpose of developing any products or services that would compete with CST products or services, (c) not alter or remove from the Products any trademarks, trade names, logos, patent or copyright notices or markings, (d) use the Products solely in accordance with CST Product Terms of Sale and any applicable documentation, and (e) comply with any license, terms of service or similar agreement with respect to any third party products or services used by Customer in connection with the Products.