GeneticistAssistant

Registration Guide

SoftGenetics, LLC ·100 Oakwood Ave, Suite 350, State College, PA 16803 · 1-888-791-1270 · www.softgenetics.com

Install the License Server Manager

- Download the License Server Manager installation file: ftp://tsupport:supportsoft@softgenetics.com/LSM/Srv%20Setup.EXE
- Right-click on the downloaded Srv%20Setup.EXE file and choose Run as Administrator to install
 - Once installation is completed, choose to Launch License Server

Register the LSM

- Open License Server Manager from the Icon Tray by clicking on LSM icon in System Tray
- Initially, you will see the message below, click OK to proceed

X

This program has not been registered. To complete the registration process, click on Help and select Register.

SoftGenetics License Server



In the LSM window that opens, go to Help > Register

Register the LSM

• Select "Geneticist Assistant" as Register Product Name.

🔌 Register Pro	duct	
Register Produc	t Name Geneticist Assistant 💌	
Register Online	Offline Registration	
Request ID	UU55RVRNeEFVTUVkak5EUKVPRIJgT	XdBREp6TWtRMFV6TTJ
Account		
Password		
Email		
D		
Hemove Lic	ense Re	egister <u>C</u> ancel

 Input your Account, Password and Email and click Register

Register the LSM

 You will be prompted with a message to restart License Server Manager for the changes to take effect, followed by a successful registration message

Message

The SoftGenetics License Server Manager needs to be restarted for changes to take affect. Please make sure that no clients are connected before restarting.



Restart Later

Information



You have register your application successfully.

X



Offline Registration

 If the LSM is unable to connect to the server online, you will receive the message below, indicating the offline registration must be used



Click OK

Offline Registration

 To proceed with offline registration, copy and paste the Request ID from the "Register Online" tab into an email and send to tech_support@softgenetics.com

🔌 Register Product 📃 💷	
Register Product Name Geneticist Assistant 💌	
Register Online Offline Registration	
Request ID EEFVTUVkak5EUkVPRJgTXdBREp6TWtRMFV6TTJjRE0yTU	Undo
Account	Cut
Password	Сору
Email	Paste
	Delete
	Select All
Remove License <u>Cance</u>	Right to left Reading order Show Unicode control characters

Offline Registration

- You will receive an email response with a registration ID
- Copy and paste the Registration ID into the field on the "Offline Registration" tab and click Register.

📚 Register Product	- • •
Register Product Name Geneticist Assistant 💌	
Register Online Offline Registration	
Register ID	Register

Configure Geneticist Assistant

 The first time you log into GA, you will be prompted to configure the license settings





Configure Geneticist Assistant

The License Setting window opens.

A Settings		1.2.6
License Settings		
Server(Hane or IP Address):		
Parts	50000	
<u></u>		
		OK Cencel

- Enter the name or IP Address where the LSM is located. Leave Port set to the default 50000
 * If it is on the same computer, enter 127.0.0.1
- Click OK



Need Assistance?

Call or Email SoftGenetics 1-888-791-1270 or 814-237-9340 tech_support@softgenetics.com