

#9998

Store at 4°C

# BSA

50 grams



**Orders** ■ 877-616-CELL (2355)  
orders@cellsignal.com

**Support** ■ 877-678-TECH (8324)  
info@cellsignal.com

**Web** ■ www.cellsignal.com

rev. 02/02/21

**For Research Use Only. Not For Use In Diagnostic Procedures.**

**Description:** Ultra pure bovine serum albumin (BSA) is used as a reagent during Western blotting. BSA can also be used in other protocols that involve antibody binding including ELISA, DELFIA, flow cytometry and immunofluorescence/immunocytochemistry.

**Storage:** Store at 4°C.

**Companion Products:**

Tris Buffered Saline with Tween 20 (TBST - 10X) #9997

Nonfat Dry Milk #9999

Phototope®-HRP Western Blot Detection System, Anti-rabbit IgG, HRP-linked Antibody #7071

Phototope®-HRP Western Blot Detection System, Anti-mouse IgG, HRP-linked Antibody #7072

Anti-rabbit IgG, HRP-linked Antibody #7074

Anti-biotin, HRP-linked Antibody #7075

Anti-mouse IgG, HRP-linked Antibody #7076

Anti-rat IgG, HRP-linked Antibody #7077

Prestained Protein Marker, Broad Range (Premixed Format) #7720

Biotinylated Protein Ladder Detection Pack #7727

20X LumiGLO® Reagent and 20X Peroxide #7003