

#13387 Store at -20°C

**TRIM33 (E1N2Z) Rabbit mAb****Cell Signaling**  
TECHNOLOGY®**Orders:** 877-616-CELL (2355)  
orders@cellsignal.com**Support:** 877-678-TECH (8324)**Web:** info@cellsignal.com  
cellsignal.com

3 Trask Lane | Danvers | Massachusetts | 01923 | USA

**For Research Use Only. Not for Use in Diagnostic Procedures.**

Applications:	Reactivity:	Sensitivity:	MW (kDa):	Source/Isotype:	UniProt ID:	Entrez-Gene Id:
WB	H	Endogenous	150	Rabbit IgG	#Q9UPN9	51592

<b>Product Usage Information</b>	<b>Application</b> Western Blotting	<b>Dilution</b> 1:1000
<b>Storage</b>	Supplied in 10 mM sodium HEPES (pH 7.5), 150 mM NaCl, 100 µg/ml BSA, 50% glycerol and less than 0.02% sodium azide. Store at -20°C. Do not aliquot the antibody.	
<b>Specificity / Sensitivity</b>	TRIM33 (E1N2Z) Rabbit mAb recognizes endogenous levels of total TRIM33 protein. Based upon sequence alignment, this antibody is predicted to react with TRIM33 isoforms A and B, but not with other TIF family members.	
<b>Species predicted to react based on 100% sequence homology:</b>	Bovine, Dog, Horse	
<b>Source / Purification</b>	Monoclonal antibody is produced by immunizing animals with a synthetic peptide corresponding to residues near the carboxy terminus of human TRIM33 protein.	
<b>Background</b>	<p>The transcriptional intermediary factor 1 (TIF1) family represents a group of proteins with multiple histone-binding domains. In humans, this family comprises four proteins, TIF1α/TRIM24, TIF1β/TRIM28/KAP1, TIF1γ/TRIM33/Ectoderm, and TIF1δ/TRIM66, which are characterized by an amino-terminal tripartite motif (TRIM) domain consisting of a RING domain, two B boxes, a coiled-coil domain, and a carboxy-terminal PHD finger and bromodomain (1). Despite their similar overall structure, these proteins have diverse roles in transcriptional regulation. TIF1α functions as a ligand-dependent nuclear receptor coregulator and more recently has been implicated in regulating p53 stability (2). TIF1β is an intrinsic component of the N-CoR1 corepressor complex and the NuRD nucleosome-remodeling complex (3) and functions as a corepressor for Kruppel-associated box (KRAB) zinc-finger transcription factors (4). Furthermore, TIF1β promotes heterochromatin-mediated gene silencing formation by serving as a cofactor for heterochromatin protein HP1 (5). TIF1δ expression is restricted to the testis and has been shown to interact with HP1γ (6).</p> <p>In contrast, the ubiquitous nuclear protein TRIM33 does not interact with either HP1 family members or chromatin-remodeling/modifying complexes. Rather, TRIM33 plays a pivotal role in signaling cascades driven by the TGF-β superfamily of ligands (7-9). A research study suggests that TRIM33 and Smad4 compete for binding to receptor phosphorylated Smad2/3 and that TRIM33-Smad2/3 and Smad4-Smad2/3 complexes complement one another in the TGF-β-dependent control of hematopoietic cell fate (9). Other studies, however, demonstrate that TRIM33 functions to repress signal relay by the TGF-β superfamily (7-8,10). Indeed, knockout of murine <i>Trim33</i> results in embryonic lethality due to upregulated Nodal signaling (10). Mechanistically, TRIM33 functions as an E3-ubiquitin ligase and promotes monoubiquitination of Smad4, a modification that impairs its ability to associate with phospho-Smad2 (8). This negative regulatory mechanism is further substantiated by the discovery that TRIM33 disrupts transcriptionally competent Smad complexes on the promoter/enhancer regions of TGF-β-responsive genes by associating with specific epigenetic marks on histone H3, which is a requirement for activating TRIM33's monoubiquitin ligase activity toward Smad4 (11). In line with the ability of TRIM33 to regulate the development of different blood cell lineages, it was shown that loss of TRIM33 expression due to epigenetic silencing of its promoter contributes to the pathogenesis of chronic myelomonocytic leukemia (12).</p>	

**Background References**

1. Meroni, G. and Diez-Roux, G. (2005) *Bioessays* 27, 1147-57.
2. Jain, A.K. and Barton, M.C. (2009) *Cell Cycle* 8, 3668-74.
3. Underhill, C. et al. (2000) *J Biol Chem* 275, 40463-70.
4. Schultz, D.C. et al. (2001) *Genes Dev* 15, 428-43.
5. Groner, A.C. et al. (2010) *PLoS Genet* 6, e1000869.
6. Khetchoumian, K. et al. (2004) *J Biol Chem* 279, 48329-41.

7. Dupont, S. et al. (2005) *Cell* 121, 87-99.
8. Dupont, S. et al. (2009) *Cell* 136, 123-35.
9. He, W. et al. (2006) *Cell* 125, 929-41.
10. Morsut, L. et al. (2010) *Development* 137, 2571-8.
11. Agricola, E. et al. (2011) *Mol Cell* 43, 85-96.
12. Aucagne, R. et al. (2011) *J Clin Invest* 121, 2361-70.

**Species Reactivity**

Species reactivity is determined by testing in at least one approved application (e.g., western blot).

**Western Blot Buffer**

**IMPORTANT:** For western blots, incubate membrane with diluted primary antibody in 5% w/v BSA, 1X TBS, 0.1% Tween® 20 at 4°C with gentle shaking, overnight.

**Applications Key**

**WB:** Western Blotting

**Cross-Reactivity Key**

**H:** human **M:** mouse **R:** rat **Hm:** hamster **Mk:** monkey **Vir:** virus **Mi:** mink **C:** chicken **Dm:** D. melanogaster  
**X:** Xenopus **Z:** zebrafish **B:** bovine **Dg:** dog **Pg:** pig **Sc:** S. cerevisiae **Ce:** C. elegans **Hr:** horse  
**GP:** Guinea Pig **Rab:** rabbit **All:** all species expected

**Trademarks and Patents**

Cell Signaling Technology is a trademark of Cell Signaling Technology, Inc.  
 SignalSilence is a registered trademark of Cell Signaling Technology, Inc.  
 XP is a registered trademark of Cell Signaling Technology, Inc.  
 All other trademarks are the property of their respective owners. Visit [cellsignal.com/trademarks](http://cellsignal.com/trademarks) for more information.

**Limited Uses**

Except as otherwise expressly agreed in a writing signed by a legally authorized representative of CST, the following terms apply to Products provided by CST, its affiliates or its distributors. Any Customer's terms and conditions that are in addition to, or different from, those contained herein, unless separately accepted in writing by a legally authorized representative of CST, are rejected and are of no force or effect.

Products are labeled with For Research Use Only or a similar labeling statement and have not been approved, cleared, or licensed by the FDA or other regulatory foreign or domestic entity, for any purpose. Customer shall not use any Product for any diagnostic or therapeutic purpose, or otherwise in any manner that conflicts with its labeling statement. Products sold or licensed by CST are provided for Customer as the end-user and solely for research and development uses. Any use of Product for diagnostic, prophylactic or therapeutic purposes, or any purchase of Product for resale (alone or as a component) or other commercial purpose, requires a separate license from CST. Customer shall (a) not sell, license, loan, donate or otherwise transfer or make available any Product to any third party, whether alone or in combination with other materials, or use the Products to manufacture any commercial products, (b) not copy, modify, reverse engineer, decompile, disassemble or otherwise attempt to discover the underlying structure or technology of the Products, or use the Products for the purpose of developing any products or services that would compete with CST products or services, (c) not alter or remove from the Products any trademarks, trade names, logos, patent or copyright notices or markings, (d) use the Products solely in accordance with CST Product Terms of Sale and any applicable documentation, and (e) comply with any license, terms of service or similar agreement with respect to any third party products or services used by Customer in connection with the Products.