

#14231 Store at 4°C

# SimpleChIP® Chromatin IP Buffers

www.cellsignal.com

1 Kit (30 immunoprecipitations)

**Support:** 877-678-TECH (8324)  
www.cellsignal.com/support

**Orders:** 877-616-CELL (2355)  
orders@cellsignal.com

rev. 06/05/17

**For Research Use Only. Not For Use In Diagnostic Procedures.**

Products Included	Product #	Quantity
ChIP Buffer (10X)	7008	20 ml
ChIP Elution Buffer (2X)	7009	7 ml
5M NaCl	7010	3 ml

**Description:** This product is offered to conveniently provide additional ChIP buffer reagents for preparing, immunoprecipitating, washing, and eluting chromatin using our SimpleChIP® (#9002, #9003) and SimpleChIP® Plus (#9004, #9005) Enzymatic Chromatin IP Kits, as well as our SimpleChIP® Plus Sonication Chromatin IP Kit (#56383). These SimpleChIP® kits provide all the reagents required for performing the recommended quantity of chromatin preparations (or optimizations) and chromatin immunoprecipitation (ChIP) assays for each kit, however there are instances where extra ChIP Buffer, ChIP Elution buffer, and NaCl are desired.

**Storage:** Store all three buffers at 4°C. They are stable for 12 months upon receipt when stored properly.

**Directions for Use:** Use as directed in the respective SimpleChIP® Chromatin IP Kit protocol.

**For product specific protocols and a complete listing of recommended companion products please see the product web page at [www.cellsignal.com](http://www.cellsignal.com)**

Thank you for your recent purchase. If you would like to provide a review visit [www.cellsignal.com/comments](http://www.cellsignal.com/comments).

© 2014 Cell Signaling Technology, Inc.  
SimpleChIP and Cell Signaling Technology are trademarks of Cell Signaling Technology, Inc.



**Applications:** W—Western IP—Immunoprecipitation IHC—Immunohistochemistry ChIP—Chromatin Immunoprecipitation IF—Immunofluorescence F—Flow cytometry E-P—ELISA-Peptide **Species Cross-Reactivity:** H—human M—mouse R—rat Hm—hamster Mk—monkey Mi—mink C—chicken Dm—D. melanogaster X—Xenopus Z—zebrafish B—bovine Dg—dog Pg—pig Sc—S. cerevisiae Ce—C. elegans Hr—Horse All—all species expected Species enclosed in parentheses are predicted to react based on 100% homology.