

ERGO Plate Stands



ERGONOMICAL PIPETTING

NEW

Dispense with maximum comfort



Ergo Plate Stands from 4titude™ are now available

- Designed for holding microplates or skirted PCR plates
- Height and angle fully adjustable
- Less stress and strain on the upper body
- Aspirate and dispense in maximum comfort
- Easy to adjust
- Also holds 1.5ml and 2.0ml microcentrifuge tubes
- Available in four colours
- Alphanumeric index (24well, 96well and 384well) cards included
- Sturdy and robust design



Available in blue, red, green and yellow for colour coordination of samples or the laboratory

Plesmanlaan 1d
2333 BZ Leiden
The Netherlands
T. +31 (0)71 568 10 00
T. Belgium: 0800 71640
F. +31 (0)71 568 10 10
info@bioke.com
www.bioke.com

BIOKÉ
sharing knowledge

ERGO Plate Stands

ERGONOMICAL PIPETTING

The height and angle of the Ergo Stand can be adjusted for maximum comfort during pipetting sequences. This reduces the stress and strain on the upper body and can help to prevent repetitive strain disorders (RSI) and work-related upper-body limb disorders (WRULD).

The stand includes eight spaces for 1.5ml or 2.0ml microcentrifuge tubes. This is especially useful for when pipetting master mixes into PCR plates.

For colour coding or identification of your work, we offer the Ergo Stands in 4 colours: red, green, blue or yellow.



Positions for 1.5ml or 2.0ml microcentrifuge tubes



Less stress and strain on the upper body



Fully adjustable height and angle for maximum pipetting comfort

Ordering Information

Code	Description	Stands/Case
4ti-0393/B	Ergo Plate Stand, Blue	1
4ti-0393/R	Ergo Plate Stand, Red	1
4ti-0393/G	Ergo Plate Stand, Green	1
4ti-0393/Y	Ergo Plate Stand, Yellow	1