

Store at -20C  
#5638

## Phospho-TNK1 (Tyr277) (D46E7) Rabbit mAb



**Cell Signaling**  
TECHNOLOGY®

**Orders:** 877-616-CELL (2355)  
orders@cellsignal.com

**Support:** 877-678-TECH (8324)

**Web:** info@cellsignal.com  
cellsignal.com

3 Trask Lane | Danvers | Massachusetts | 01923 | USA

**For Research Use Only. Not for Use in Diagnostic Procedures.**

Applications:	Reactivity:	Sensitivity:	MW (kDa):	Source/Isotype:	UniProt ID:	Entrez-Gene Id:
WB	H	Endogenous	72 TNK1, 58 TNK1-C17orf61	Rabbit IgG	#Q13470	8711

### Product Usage Information

#### Application

Western Blotting

#### Dilution

1:1000

### Storage

Supplied in 10 mM sodium HEPES (pH 7.5), 150 mM NaCl, 100 µg/ml BSA, 50% glycerol and less than 0.02% sodium azide. Store at -20°C. Do not aliquot the antibody.

### Specificity / Sensitivity

Phospho-TNK1 (Tyr277) (D46E7) Rabbit mAb detects endogenous levels of TNK1 protein only when the amplified gene product is phosphorylated at Tyr277.

### Species predicted to react based on 100% sequence homology:

Mouse, Rat

### Source / Purification

Monoclonal antibody is produced by immunizing animals with a synthetic phosphopeptide corresponding to residues surrounding Tyr277 of human TNK1 protein.

### Background

Tyrosine kinase non-receptor 1 (TNK1) is related to the Ack1 (TNK2) non-receptor kinase that binds cdc42 and inhibits GTPase activity of this cell cycle regulator. TNK1 is broadly expressed in embryogenic tissues and leukemia cell lines, but is restricted to select adult tissues (1). TNK1 is a putative 72 kDa protein comprised of an N-terminal kinase domain, a central SH3 domain and a proline-rich tail. Interaction with PLCγ *in vitro* indicates a possible role in phospholipid signal transduction pathways (2). Though the exact mechanism is currently unclear, active TNK1 may play a role in regulating cell death by preventing TNF-α-induced NF-κB transactivation (3). Phosphorylation of TNK1 on Tyr277 was identified at Cell Signaling Technology (CST) using PhosphoScan®, CST's LC-MS/MS platform for phosphorylation site discovery (4) and also reported independently in another publication using MS technology (5). Phosphorylation of TNK1 at Tyr277 was observed in select carcinoma cell lines and in tumors. A constitutively active, 58 kDa truncated TNK1 kinase resulting from fusion between the TNK1 and C17orf61 genes is seen in some cells (5). For additional information, visit PhosphoSitePlus™, CST's modification site knowledgebase, at [www.phosphosite.org](http://www.phosphosite.org).

### Background References

1. Hoehn, G.T. et al. (1996) *Oncogene* 12, 903-13.
2. Felschow, D.M. et al. (2000) *Biochem Biophys Res Commun* 273, 294-301.
3. Azoitei, N. et al. (2007) *Oncogene* 26, 6536-45.
4. Rush, J. et al. (2005) *Nat Biotechnol* 23, 94-101.
5. Gu, T.L. et al. (2010) *Leukemia* 24, 861-5.

### Species Reactivity

Species reactivity is determined by testing in at least one approved application (e.g., western blot).

### Western Blot Buffer

IMPORTANT: For western blots, incubate membrane with diluted primary antibody in 5% w/v BSA, 1X TBS, 0.1% Tween® 20 at 4°C with gentle shaking, overnight.

### Applications Key

**WB:** Western Blotting

### Cross-Reactivity Key

**H:** human **M:** mouse **R:** rat **Hm:** hamster **Mk:** monkey **Vir:** virus **Mi:** mink **C:** chicken **Dm:** D. melanogaster  
**X:** Xenopus **Z:** zebrafish **B:** bovine **Dg:** dog **Pg:** pig **Sc:** S. cerevisiae **Ce:** C. elegans **Hr:** horse  
**GP:** Guinea Pig **Rab:** rabbit **All:** all species expected

**Trademarks and Patents**

Cell Signaling Technology is a trademark of Cell Signaling Technology, Inc.

All other trademarks are the property of their respective owners. Visit [cellsignal.com/trademarks](https://cellsignal.com/trademarks) for more information.

**Limited Uses**

Except as otherwise expressly agreed in a writing signed by a legally authorized representative of CST, the following terms apply to Products provided by CST, its affiliates or its distributors. Any Customer's terms and conditions that are in addition to, or different from, those contained herein, unless separately accepted in writing by a legally authorized representative of CST, are rejected and are of no force or effect.

Products are labeled with For Research Use Only or a similar labeling statement and have not been approved, cleared, or licensed by the FDA or other regulatory foreign or domestic entity, for any purpose. Customer shall not use any Product for any diagnostic or therapeutic purpose, or otherwise in any manner that conflicts with its labeling statement. Products sold or licensed by CST are provided for Customer as the end-user and solely for research and development uses. Any use of Product for diagnostic, prophylactic or therapeutic purposes, or any purchase of Product for resale (alone or as a component) or other commercial purpose, requires a separate license from CST. Customer shall (a) not sell, license, loan, donate or otherwise transfer or make available any Product to any third party, whether alone or in combination with other materials, or use the Products to manufacture any commercial products, (b) not copy, modify, reverse engineer, decompile, disassemble or otherwise attempt to discover the underlying structure or technology of the Products, or use the Products for the purpose of developing any products or services that would compete with CST products or services, (c) not alter or remove from the Products any trademarks, trade names, logos, patent or copyright notices or markings, (d) use the Products solely in accordance with CST Product Terms of Sale and any applicable documentation, and (e) comply with any license, terms of service or similar agreement with respect to any third party products or services used by Customer in connection with the Products.