

#8595 Store at -20°C

## Nuclear Receptor Antibody Sampler Kit

1 Kit (8 x 20 microliters)



**Cell Signaling**  
TECHNOLOGY®

**Orders:** 877-616-CELL (2355)  
orders@cellsignal.com

**Support:** 877-678-TECH (8324)

**Web:** info@cellsignal.com  
cellsignal.com

3 Trask Lane | Danvers | Massachusetts | 01923 | USA

**For Research Use Only. Not for Use in Diagnostic Procedures.**

Product Includes	Product #	Quantity	Mol. Wt	Isotype/Source
RAR $\alpha$ (E6Z6K) Rabbit mAb	62294	20 $\mu$ l	60 kDa	Rabbit IgG
RAR $\gamma$ 1 (D3A4) XP® Rabbit mAb	8965	20 $\mu$ l	58 kDa	Rabbit IgG
RXR $\alpha$ (D6H10) Rabbit mAb	3085	20 $\mu$ l	53 kDa	Rabbit IgG
Glucocorticoid Receptor (D8H2) XP® Rabbit mAb	3660	20 $\mu$ l	80, 91, 94 kDa	Rabbit IgG
Progesterone Receptor A/B (D8Q2J) XP® Rabbit mAb	8757	20 $\mu$ l	90 (PR-A), 118 (PR-B) kDa	Rabbit IgG
Androgen Receptor (D6F11) XP® Rabbit mAb	5153	20 $\mu$ l	110 kDa	Rabbit IgG
Estrogen Receptor $\alpha$ (D8H8) Rabbit mAb	8644	20 $\mu$ l	66 kDa	Rabbit IgG
PPAR $\gamma$ (C26H12) Rabbit mAb	2435	20 $\mu$ l	53, 57 kDa	Rabbit IgG
Anti-rabbit IgG, HRP-linked Antibody	7074	100 $\mu$ l		Goat

Please visit [cellsignal.com](https://www.cellsignal.com) for individual component applications, species cross-reactivity, dilutions, protocols, and additional product information.

### Description

The Nuclear Receptor Antibody Sampler Kit provides an economical means to evaluate the presence and status of nuclear receptors. This kit contains enough primary antibody to perform two western blots per primary.

### Storage

Supplied in 10 mM sodium HEPES (pH 7.5), 150 mM NaCl, 100  $\mu$ g/ml BSA, 50% glycerol and less than 0.02% sodium azide. Store at -20°C. Do not aliquot the antibody.

### Background

Nuclear Receptors are transcription factors responsible for sensing bioactive molecules, including steroid and thyroid hormones. They are regulated by multiple posttranslational modifications, which in turn impacts their ability to regulate the expression of specific genes involved in the control of reproduction, development, and metabolism.

Androgen receptor (AR), a zinc finger transcription factor belonging to the nuclear receptor superfamily, is activated by phosphorylation and dimerization upon ligand binding (1). This promotes nuclear localization and binding of AR to androgen response elements in androgen target genes. AR plays a crucial role in several stages of male development and the progression of prostate cancer (2,3).

Estrogen receptor  $\alpha$  (ER $\alpha$ ), a member of the steroid receptor superfamily, contains highly conserved DNA binding and ligand binding domains (4). Through its estrogen-independent and estrogen-dependent activation domains (AF-1 and AF-2, respectively), ER $\alpha$  regulates transcription by recruiting coactivator proteins and interacting with general transcriptional machinery (5).

Glucocorticoid hormones control cellular proliferation, inflammation, and metabolism through their association with the glucocorticoid receptor (GR)/NR3C1, a member of the nuclear hormone receptor superfamily of transcription factors (6).

Peroxisome proliferator-activated receptor  $\gamma$  (PPAR $\gamma$ ) is a member of the ligand-activated nuclear receptor superfamily and functions as a transcriptional activator (7). PPAR $\gamma$  is preferentially expressed in adipocytes, as well as in vascular smooth muscle cells and macrophages (8). Besides its role in mediating adipogenesis and lipid metabolism (8), PPAR $\gamma$  also modulates insulin sensitivity, cell proliferation, and inflammation (9).

Human progesterone receptor (PR) is expressed as two forms: the full length PR B and the short form PR A. PR A lacks the first 164 amino acid residues of PR B (10,11). Both PR A and PR B are ligand activated, but differ in their relative ability to activate target gene transcription (12,13).

Nuclear retinoic acid receptors (RARs) consist of three subtypes encoded by separate genes:  $\alpha$  (NR1B1),  $\beta$  (NR1B2), and  $\gamma$  (NR1B3). For each subtype, there are at least two isoforms, which are generated by differential promoter usage and alternative splicing and differ only in their N-terminal regions. Retinoids,

which are metabolites of vitamin A, serve as ligands for RARs (14). RARs function as ligand-dependent transcriptional regulators and are found to be heterodimerized with retinoid X receptors (RXRs). These transcriptionally active dimers regulate the expression of genes involved in cellular differentiation, proliferation, and apoptosis (15,16).

The human retinoid X receptors are encoded by three distinct genes (*RXRα*, *RXRβ*, and *RXRγ*) and bind selectively and with high affinity to the vitamin A derivative, 9-*cis*-retinoic acid. RXRs are type-II nuclear hormone receptors that are largely localized to the nuclear compartment independent of ligand binding. Nuclear RXRs form heterodimers with nuclear hormone receptor subfamily 1 proteins, including thyroid hormone receptor, retinoic acid receptors, vitamin D receptor, peroxisome proliferator-activated receptors, liver X receptors, and farnesoid X receptor (17).

## Background References

1. Li, J. and Al-Azzawi, F. (2009) *Maturitas* 63, 142-8.
2. Avila, D.M. et al. (2010) *J Steroid Biochem Mol Biol* 76, 135-42.
3. Montgomery, J.S. et al. (2001) *J Pathol* 195, 138-46.
4. Mangelsdorf, D.J. et al. (1995) *Cell* 83, 835-9.
5. Glass, C.K. and Rosenfeld, M.G. (2000) *Genes Dev* 14, 121-41.
6. Yamamoto, K.R. (1985) *Annu Rev Genet* 19, 209-52.
7. Tontonoz, P. et al. (1995) *Curr Opin Genet Dev* 5, 571-6.
8. Rosen, E.D. et al. (1999) *Mol Cell* 4, 611-7.
9. Murphy, G.J. and Holder, J.C. (2000) *Trends Pharmacol Sci* 21, 469-74.
10. Evans, R.M. (1988) *Science* 240, 889-95.
11. Kastner, P. et al. (1990) *EMBO J* 9, 1603-14.
12. Giangrande, P.H. et al. (2000) *Mol Cell Biol* 20, 3102-15.
13. Wen, D.X. et al. (1994) *Mol Cell Biol* 14, 8356-64.
14. Rochette-Egly, C. and Germain, P. (2009) *Nucl Recept Signal* 7, e005.
15. Delacroix, L. et al. (2010) *Mol Cell Biol* 30, 231-44.
16. Eifert, C. et al. (2006) *Mol Reprod Dev* 73, 796-824.
17. Gronemeyer, H. et al. (2004) *Nat Rev Drug Discov* 3, 950-64.

## Trademarks and Patents

Cell Signaling Technology is a trademark of Cell Signaling Technology, Inc. U.S. Patent No. 7,429,487, foreign equivalents, and child patents deriving therefrom. All other trademarks are the property of their respective owners. Visit [cellsignal.com/trademarks](http://cellsignal.com/trademarks) for more information.

## Limited Uses

Except as otherwise expressly agreed in a writing signed by a legally authorized representative of CST, the following terms apply to Products provided by CST, its affiliates or its distributors. Any Customer's terms and conditions that are in addition to, or different from, those contained herein, unless separately accepted in writing by a legally authorized representative of CST, are rejected and are of no force or effect.

Products are labeled with For Research Use Only or a similar labeling statement and have not been approved, cleared, or licensed by the FDA or other regulatory foreign or domestic entity, for any purpose. Customer shall not use any Product for any diagnostic or therapeutic purpose, or otherwise in any manner that conflicts with its labeling statement. Products sold or licensed by CST are provided for Customer as the end-user and solely for research and development uses. Any use of Product for diagnostic, prophylactic or therapeutic purposes, or any purchase of Product for resale (alone or as a component) or other commercial purpose, requires a separate license from CST. Customer shall (a) not sell, license, loan, donate or otherwise transfer or make available any Product to any third party, whether alone or in combination with other materials, or use the Products to manufacture any commercial products, (b) not copy, modify, reverse engineer, decompile, disassemble or otherwise attempt to discover the underlying structure or technology of the Products, or use the Products for the purpose of developing any products or services that would compete with CST products or services, (c) not alter or remove from the Products any trademarks, trade names, logos, patent or copyright notices or markings, (d) use the Products solely in accordance with CST Product Terms of Sale and any applicable documentation, and (e) comply with any license, terms of service or similar agreement with respect to any third party products or services used by Customer in connection with the Products.