

#9920 Store at -20°C

## Phospho-(Ser/Thr) Kinase Substrate Antibody Sampler Kit

1 Kit (6 x 20 microliters)



**Cell Signaling**  
TECHNOLOGY®

**Orders:** 877-616-CELL (2355)  
orders@cellsignal.com

**Support:** 877-678-TECH (8324)

**Web:** info@cellsignal.com  
cellsignal.com

3 Trask Lane | Danvers | Massachusetts | 01923 | USA

**For Research Use Only. Not for Use in Diagnostic Procedures.**

Product Includes	Product #	Quantity	Mol. Wt	Isotype/Source
Phospho-AMPK Substrate Motif [LXRX(pS/pT)] MultiMab® Rabbit mAb mix	5759	20 µl		Rabbit
Phospho-Akt Substrate (RXXS*/T*) (110B7E) Rabbit mAb	9614	20 µl		Rabbit IgG
Phospho-PKA Substrate (RRXS*/T*) (100G7E) Rabbit mAb	9624	20 µl		Rabbit IgG
Phospho-ATM/ATR Substrate Motif [(pS/pT) QG] MultiMab® Rabbit mAb mix	6966	20 µl		Rabbit IgG
Phospho-PKC Substrate Motif [(R/K)XpSX(R/K)] MultiMab® Rabbit mAb mix	6967	20 µl		Rabbit IgG
Phospho-CDK Substrate Motif [(K/H)pSP] MultiMab® Rabbit mAb mix	9477	20 µl		Rabbit IgG
Anti-rabbit IgG, HRP-linked Antibody	7074	100 µl		Goat

Please visit [cellsignal.com](http://cellsignal.com) for individual component applications, species cross-reactivity, dilutions, protocols, and additional product information.

### Description

The Phospho-(Ser/Thr) Kinase Substrate Antibody Sampler Kit provides an economical means to investigate the downstream activity of select serine/threonine kinases. The kit contains enough primary antibody to perform two western blot experiments with each antibody.

### Storage

Supplied in 10 mM sodium HEPES (pH 7.5), 150 mM NaCl, 100 µg/ml BSA, 50% glycerol and less than 0.02% sodium azide. Store at -20°C. Do not aliquot the antibody.

### Background

Phospho-(Ser/Thr) kinases and phosphatases play critical roles in a wide range of biological processes. Each phospho-(Ser/Thr) kinase phosphorylates serine or threonine within a specific motif. Akt phosphorylates substrates at a serine or threonine only in a conserved motif characterized by arginine at positions -5 and -3 (1). Conventional PKC isozymes phosphorylate substrates containing serine or threonine, with arginine or lysine at the -3, -2 and +2 positions, and a hydrophobic amino acid at position +1 (2,3). A consensus phosphorylation site of PKA is serine or threonine with arginine at the -2 and -3 positions (3). AMPK phosphorylates consensus motif (L/M)XRX(S/T)XXXL (6). Antibodies recognizing the LXRX(S/T) motif are very useful in the identification of AMPK substrates. The consensus amino acid sequence for CDK substrate is (K/R)(S\*)PX(K/R), where denotes any one of the 20 amino acids and S\* is the phosphorylation site (4-6). ATM and the related kinase ATR phosphorylate serine or threonine in an S\*/T\*Q motif (7,8).

Antibodies specific to particular kinase substrates are invaluable reagents in determining kinase activity and identifying potential new kinase substrates. CST has developed antibodies that recognize phosphorylated serine or threonine within the context of a protein motif that is phosphorylated by Akt, PKC, PKA, MAPK/CDK, CDKs or ATM/ATR. As shown by peptide pairing ELISA, each phospho-(Ser/Thr) kinase substrate antibody in this sampler kit is specific to its kinase substrate motif.

### Background References

1. Alessi, D.R. et al. (1996) *FEBS Lett* 399, 333-8.
2. Nishikawa, K. et al. (1997) *J Biol Chem* 272, 952-60.
3. Pearson, R.B. and Kemp, B.E. (1991) *Methods Enzymol* 200, 62-81.
4. Songyang, Z. et al. (1996) *Mol Cell Biol* 16, 6486-93.
5. Songyang, Z. (1999) *Prog Biophys Mol Biol* 71, 359-72.
6. Holmes, J.K. and Solomon, M.J. (1996) *J Biol Chem* 271, 25240-6.
7. Kastan, M.B. and Lim, D.S. (2000) *Nat Rev Mol Cell Biol* 1, 179-86.
8. Zhao, H. and Piwnicka-Worms, H. (2001) *Mol Cell Biol* 21, 4129-39.

### Trademarks and Patents

Cell Signaling Technology is a trademark of Cell Signaling Technology, Inc. All other trademarks are the property of their respective owners. Visit [cellsignal.com/trademarks](http://cellsignal.com/trademarks) for more information.

### Limited Uses

Except as otherwise expressly agreed in a writing signed by a legally authorized representative of CST, the following terms apply to Products provided by CST, its affiliates or its distributors. Any Customer's terms and conditions that are in addition to, or different from, those contained herein, unless separately accepted in

writing by a legally authorized representative of CST, are rejected and are of no force or effect.

Products are labeled with For Research Use Only or a similar labeling statement and have not been approved, cleared, or licensed by the FDA or other regulatory foreign or domestic entity, for any purpose. Customer shall not use any Product for any diagnostic or therapeutic purpose, or otherwise in any manner that conflicts with its labeling statement. Products sold or licensed by CST are provided for Customer as the end-user and solely for research and development uses. Any use of Product for diagnostic, prophylactic or therapeutic purposes, or any purchase of Product for resale (alone or as a component) or other commercial purpose, requires a separate license from CST. Customer shall (a) not sell, license, loan, donate or otherwise transfer or make available any Product to any third party, whether alone or in combination with other materials, or use the Products to manufacture any commercial products, (b) not copy, modify, reverse engineer, decompile, disassemble or otherwise attempt to discover the underlying structure or technology of the Products, or use the Products for the purpose of developing any products or services that would compete with CST products or services, (c) not alter or remove from the Products any trademarks, trade names, logos, patent or copyright notices or markings, (d) use the Products solely in accordance with CST Product Terms of Sale and any applicable documentation, and (e) comply with any license, terms of service or similar agreement with respect to any third party products or services used by Customer in connection with the Products.